

## GeneMAGNET Human and Animal Tissue DNA Purification Kit

Kit for isolation of total DNA from human and animal tissues

○ **Cat. no. E3425**

EURx Ltd. 80-297 Gdansk Poland  
ul. Przyrodników 3, NIP 957-07-05-191  
KRS 0000202039, [www.eurx.com.pl](http://www.eurx.com.pl)  
orders: email: [orders@eurx.com.pl](mailto:orders@eurx.com.pl)  
tel. +48 58 524 06 97, fax +48 58 341 74 23



# Table of Contents

Introductory Notes.....	3
Equipment and reagents to be supplied by the user .....	3
Protocol.....	4
Part I Sample preparation .....	4
A. Solid tissues/Insects/Rodent tails .....	4
B. Paraffin-embedded tissues.....	4
C. Formalin-fixed tissues .....	5
D. Liquid tissues/body fluids.....	5
E. Saliva/swab .....	6
F. Semen .....	6
G. Cultured cells .....	6
H. Hair.....	7
I. Urine .....	7
Part II DNA isolation.....	8
Appendix 1 DNA isolation from bacterial cells.....	9
Safety Information .....	10

Content	96 preps E3425-01	Storage/Stability
Lyse T	30 ml	15-25°C
RNase A (10 mg/ml)	0.22 ml	2-8°C
Proteinase K (20 mg/ml)	2.2 ml	-20°C
Sol T	30 ml	15-25°C
Wash T1	2 x 88 ml	15-25°C
Wash T2	88 ml	15-25°C
Elution	24 ml	15-25°C
Magnetic Beads	40 ml	2-8°C
Protocol	1	

# Introductory Notes

**NOTE 1 • Kit Specification.** The kit is designed for isolation of total DNA (genomic, mitochondrial) from a variety of tissues and biological liquids by silica-covered Magnetic Beads. The kit is designed for manual or automatic use.

**NOTE 2 • Maximum Sample Amount.** One preparation enables purification of DNA from up to 25 mg solid tissues or 100 µl liquid tissues. The volume of the sample can be scaled up or down with the proportion 2-2.5 : 2 : 3.5 (sample with lysis bufer, Sol T and Magnetic Beads).

**NOTE 3 • Kit Compounds Storage.** Once the kit is unpacked, store components at room temperature, with the exception of Magnetic Beads, RNase A and Proteinase K. Magnetic Beads and RNase A should be kept at 2-8°C and Proteinase K at -20°C.

**NOTE 4 • Internal control.** Commercially available amplification systems may require the introduction of an internal control into the purification procedure. Internal control DNA should be added to the lysis buffer Lyse T.

**NOTE 5 • Additional recommendations.** Avoid cross contamination of samples during isolation procedure by changing pipette tips whenever necessary.

**NOTE 6 • DNA purification efficiency.** Working with Magnetic Beads, depending on the source material, the user can obtain high DNA concentrations of the isolates. PCR optimization is required in order to avoid nonspecific products or low reaction efficiency. The optimal genomic DNA template is 5-50 ng per 50 µl PCR reaction volume.

**NOTE 7 • Maintaining Good Working Practice.** All solutions should be kept tightly closed to avoid evaporation and resulting concentration changes of buffer components. To obtain high quality DNA, stick carefully to the protocol provided below.

## *Equipment and reagents to be supplied by the user*

- Magnetic stand E0361 for 16 tubes, E0362 for 24 tubes, E0363 for 96-well plate. To be purchased separately.
- Microcentrifuge, disposable gloves, sterile pipette tips, sterile 1.5-2 ml tubes or 96-well plates with working volume at least 800 µl, a heating block capable of incubation at 37-70°C.
- [1 M] Dithiothreitol (DTT), ethanol [96 -100% v/v], xylene, PBS (cat. no E0281) and 10 mM Tris pH 8.0 (add 10 ml 1M Tris-HCl pH 8.0 EURx cat. no E0273 to 990 ml sterile water). Optional AFR01 (E0328) antifoam agent for foam reduction during solid tissues homogenization with mechanical homogenizer.



# Protocol

## Part I Sample preparation

### A. Solid tissues/Insects/Rodent tails

1. Choose tissue homogenization method.
  - a) Grind tissue fragment under liquid nitrogen to a fine powder using previously cooled mortar and pestle. Place sample material (up to 25 mg) in 2 ml Eppendorf tube and centrifuge the powder to the bottom of the tube. Add 200  $\mu$ l of buffer **Lyse T** and suspend the precipitate thoroughly.
    - *To obtain high yield of DNA a tissue fragment should be thoroughly grinded to a fine powder.*
  - b) Place up to 25 mg of tissue in 2 ml Eppendorf tube. Add 200  $\mu$ l **Lyse T** and homogenize the sample using a mechanical homogenizer.
    - *Lyse T buffer contains detergents which produce foam during rapid homogenization with mechanical homogenizer. Optionally add AFR01 antifoam agent.*
    - *Alternatively homogenize the sample in 100  $\mu$ l PBS and add 100  $\mu$ l Lyse T.*
  - c) Cut tissue fragment (up to 25 mg) into small pieces. Place the sample in 2 ml Eppendorf tube and add 200  $\mu$ l of buffer **Lyse T**.
2. Add 2  $\mu$ l of **RNase A** and 20  $\mu$ l of **Proteinase K**. Mix by inverting or vortexing the tube.
3. Incubate at 56°C until the tissue is completely lysed (at least 1-3 h). Mix by inverting or vortexing every 15-30 min.
  - *Samples can be lysed overnight, if needed.*
4. Spin down with maximal speed for 3 min and transfer the supernatant into the new tube.
  - *This step enables to remove undigested tissue fragments, hair and bones.*
5. Follow the point 1. Part II. of the DNA isolation protocol.

### B. Paraffin-embedded tissues

1. Prepare a small section (up to 25 mg) from block of embedded tissue. Place the sample in 2 ml Eppendorf tube.
2. Add 1 ml xylene. Vortex vigorously. Incubate at room temperature for 15 min.
3. Centrifuge for 3 min at maximal speed. Remove supernatant by pipetting.
4. Add 1 ml xylene to the pellet, mix by vortexing.
5. Centrifuge for 3 min at maximal speed. Remove supernatant by pipetting.

6. Add 1 ml ethanol (96–100%) to the pellet. Mix by vortexing or inverting the tube.
7. Centrifuge for 3 min at maximum speed. Remove supernatant by pipetting.
8. Repeat steps 6–7 once.
9. Incubate the open tube at 37°C until the ethanol has evaporated (app. 15 min).
10. Resuspend the tissue pellet in 200 µl of buffer **Lyse T**.
11. Add 2 µl of **RNase A** and 20 µl of **Proteinase K**. Mix by inverting or vortexing the tube.
12. Incubate at 56°C until the tissue is completely lysed (at least 1-3 h). Mix by inverting or vortexing every 15-30 min.
  - *Samples can be lysed overnight, if needed.*
13. Follow the point 1. Part II. of the DNA isolation protocol.

#### C. Formalin-fixed tissues

1. Wash tissue sample twice with PBS to remove fixative. Discard PBS.
2. Cut tissue fragment (up to 25 mg) into small pieces. Place the sample in 2 ml Eppendorf tube and add 200 µl of buffer **Lyse T**.
3. Add 2 µl of **RNase A** and 20 µl of **Proteinase K**. Mix by inverting or vortexing the tube.
4. Incubate at 56°C until the tissue is completely lysed (at least 1-3 h). Mix by inverting or vortexing every 15-30 min.
  - *Samples can be lysed overnight, if needed.*
5. Follow the point 1. Part II. of the DNA isolation protocol.

#### D. Liquid tissues/body fluids

(blood, saliva, plasma, serum, brain-spinal cord liquid among others)

1. To 100 µl of sample add 2 µl **RNase A** and 100 µl **Lyse T**.
  - *For sample volumes less than 100 µl, add PBS to adjust the volume to 100 µl.*
2. Mix thoroughly by vortexing the tube.
3. Incubate for 5 min at room temperature.
4. Add 10 µl **Proteinase K**.
5. Follow the point 1. Part II. of the DNA isolation protocol.



## E. Saliva/swab

1. To 100  $\mu\text{l}$  of sample add 2  $\mu\text{l}$  **RNase A** and 100  $\mu\text{l}$  **Lyse T**.
  - In case of dry swabs: add 300  $\mu\text{l}$  PBS into the sample, mix and incubate for at least 5 min. Transfer 100  $\mu\text{l}$  of sample to the new tube.
  - In case of swabs in storage buffer (with guanidine) transfer 100  $\mu\text{l}$  of sample to the new tube.
  - It is possible to add 200  $\mu\text{l}$  of Lyse T, 2  $\mu\text{l}$  RNase A and 10  $\mu\text{l}$  Proteinase K directly to the dry swab, mix, incubate for 5 min at room temperature and follow the point 1. Part II of the protocol. Note that the dry swab can absorb some of the liquid.
2. Mix thoroughly by vortexing the tube.
3. Incubate for 5 min at room temperature.
4. Add 10  $\mu\text{l}$  **Proteinase K**.
5. Follow the point 1. Part II. of the DNA isolation protocol.

## F. Semen

1. Add 1.5 ml 10 mM Tris pH 8.0 to 100  $\mu\text{l}$  of semen, mix for 10 min by vortexing.
2. Spin down for 5 min with 10 000 x g.
3. Remove the supernatant and resuspend the pellet in 200  $\mu\text{l}$  **Lyse T**.
4. Add 2  $\mu\text{l}$  **RNase A**, 20  $\mu\text{l}$  **Proteinase K** and 10  $\mu\text{l}$  1 M DTT.
5. Mix and incubate at 56°C until the tissue is completely lysed (at least 1 h).
  - Samples can be lysed overnight, if needed.
6. Follow the point 1. Part II. of the DNA isolation protocol.

## G. Cultured cells

1. Centrifuge the cell culture (up to  $5 \times 10^6$  cells) in the 1.5–2 ml Eppendorf tube for 3 min at 1 000 x g.
2. Carefully discard the supernatant. Add to the pellet 200  $\mu\text{l}$  **Lyse T** buffer and 2  $\mu\text{l}$  **RNase A**. Suspend the cells thoroughly by vortexing for 20 sec.
3. Incubate for 5 min at room temperature.
4. Add 10  $\mu\text{l}$  **Proteinase K**.
5. Follow the point 1. Part II. of the DNA isolation protocol.

## H. Hair

1. Cut off the hair roots from the hair sample (up to 100 roots or 25 mg). Place them in the 2 ml Eppendorf tube. Add 200  $\mu$ l of buffer **Lyse T**, 20  $\mu$ l of 1M DTT and 20  $\mu$ l of **Proteinase K**. Mix by vortexing.
  - *If the hair sample doesn't contain the roots cut the hair stems into short pieces not longer than 0.5 cm.*
  - *The hair stem is the dead part of hair that contain small quantities of degraded DNA. The recommended amplicon length for PCR analysis of DNA from the hair stems is <200 bp.*
2. Incubate at 56°C until the hair sample is completely lysed (6-8 h or overnight).
3. Mix by vortexing every 1-2 h or use a shaking water bath.
4. Spin down for 3 min with the maximal speed.
  - *This step removes undigested tissue fragments.*
5. Follow the point 1. Part II. of the DNA isolation protocol.

## I. Urine

1. Add 2 ml of urine to the 2 ml Eppendorf tube.
2. Centrifuge urine in microcentrifuge for 2 min at 6 000 x g.
3. Carefully discard the supernatant without disturbing the pellet. Add to the pellet 200  $\mu$ l **Lyse T** buffer and 10  $\mu$ l **Proteinase K**.
4. Vortex for 15 sec.
5. Incubate for 60 min at 56°C, mix by inverting the tube every 15 min.
6. Follow the point 1. Part II. of the DNA isolation protocol.

## Part II DNA isolation

1. Add 200  $\mu$ l of buffer **Sol T** and mix thoroughly by vortexing or several-fold inverting.
  - o *If working on 96-well plate format transfer the samples to the wells of the plate.*
2. Resuspend **Magnetic Beads** before removing them from the storage tube by vortexing. Add 350  $\mu$ l of resuspended **Magnetic Beads** to the sample and mix by vortexing or pipetting for 1 min.
3. Separate the **Magnetic Beads** against the side of the tubes/wells. Wait until all the beads have been attached to the magnets.
4. Remove and discard the supernatant by pipetting. Remove the magnetic stand/transfer tubes to the laboratory rack, add 800  $\mu$ l of **Wash T1** and mix by pipetting or vortexing for 10 s.
5. Separate the **Magnetic Beads** against the side of the tubes/wells. Wait until all the beads have been attached to the magnets.
6. Remove and discard the supernatant by pipetting. Remove the magnetic stand/transfer tubes to the laboratory rack, add 800  $\mu$ l of **Wash T1** and mix by pipetting or vortexing for 10 s.
7. Separate the **Magnetic Beads** against the side of the tubes/wells. Wait until all the beads have been attached to the magnets.
8. Remove and discard the supernatant by pipetting. Remove the magnetic stand/transfer tubes to the laboratory rack, add 800  $\mu$ l of **Wash T2** and mix by pipetting or vortexing for 10 s.
9. Separate the **Magnetic Beads** against the side of the tubes/wells. Wait until all the beads have been attached to the magnets.
10. Remove and discard the supernatant by pipetting. Leave the open tubes/plate in magnetic stand and air dry the beads for 15-20 min.
  - o *Wash T3 contains alcohol, make sure all the solution evaporates before proceeding to the next step.*
11. Add 50-200  $\mu$ l **RNase-free water** to the tube/well and mix by pipetting or vortexing for 2 min.
12. Separate the **Magnetic Beads** against the side of the wells. After all the beads have been attached to the magnets transfer the supernatant containing the purified DNA to a suitable tube/plate. DNA is ready for analysis/manipulations. Isolated DNA can be stored either at 2-8° or at -20°.



## Appendix 1 DNA isolation from bacterial cells

**NOTE 1** • It is recommended to use this protocol for genomic DNA isolation from bacteria cultures, suspensions and tissues/body fluids containing bacteria.

**NOTE 2** • Lysozyme containing BL buffer (E0309-02) is not included in this kit, it has to be ordered separately. While ordering please indicate the buffer name and catalog number. Certain bacterial species are resistant to lysis, thus supplementary enzymes other than lysozyme may be necessary. For example, lysis of *Staphylococcus* is much more efficient with lysostaphin.

1. Mix in 1.5 ml Eppendorf tube:
  - a) 80 µl overnight bacterial culture and 100 µl **Lyse T** buffer.
  - b) Pick bacterial colony directly from Petri dish and suspend in 180 µl buffer **Lyse BG**.
  - c) Pellet bacteria from 0.1–1.5 ml overnight culture by centrifugation and discard the supernatant, ensuring that all liquid is completely removed. Resuspend the bacterial pellet in 180 µl buffer **Lyse T**.
    - For high yield isolation it is critical to completely resuspend bacterial cells.
    - The highest quality DNA is obtained from bacterial culture, which are either in log phase or early stationary phase.
2. Add 50 µl buffer **BL** and 2 µl **RNase A** to the suspension cell (point 1). Mix by several-fold inverting or vortex 3 sec.
  - For efficient lysis of some bacterial species, enzymes other than lysozyme may be necessary. Use the appropriate enzyme (with buffer BL) for the particular species.
3. Incubate the sample at 37°C for 15 min.
4. Add 20 µl **Proteinase K** to the resuspended cell pellet. Mix by several-fold inverting or vortex 3 sec.
5. Incubate the sample at 56°C for 30 min.
6. Follow the point 1. Part II. of the DNA isolation protocol.

# Safety Information

## Lyse T

---

### Warning



**H319** Causes serious eye irritation.

**P280** Wear protective gloves/protective clothing/eye protection/face protection.

**P305+P351+P338** If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**P337+P313** If eye irritation persists: Get medical advice/ attention.

## Proteinase K

---

### Danger



**H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled.

**P261** Avoid breathing vapours/spray.

**P304+P340** If inhaled: remove person to fresh air and keep comfortable for breathing.

**P342+P311** If experiencing respiratory symptoms: call a poison center or doctor/ physician.

## Sol T

---

### Warning



**H302+H332** Harmful if swallowed or if inhaled.

**H315** Causes skin irritation.

**H319** Causes serious eye irritation.

**P261** Avoid breathing vapours/spray.

**P280** Wear protective gloves/protective clothing/eye protection/face protection.

**P301+P312** If swallowed: call a poison center/ doctor if you feel unwell.

**P304+P340** If inhaled: remove person to fresh air and keep comfortable for breathing.

**P305+P351+P338** If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**P333+P313** If skin irritation or rash occurs: get medical advice/attention.

**P337+P313** If eye irritation persists: get medical advice/ attention.

**EUH208** Contains ethylenediammonium dichloride. May produce an allergic reaction.

## Wash T1 / Wash T2

---

### Danger



**H225** Highly flammable liquid and vapour.

**H319** Causes serious eye irritation.

**P210** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.



**P280** Wear protective gloves/protective clothing/eye protection/face protection.

**P305+P351+P338** If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**P403+P235** Store in a well-ventilated place. Keep cool.

**P337+P313** If eye irritation persists: get medical advice/ attention.

- **GeneMAGNET Human and Animal Tissue DNA Purification Kit is designed for rapid manual or automatic purification of total DNA (genomic, mitochondrial) from a variety of tissues and biological liquids. Purified DNA is free of contaminants, such as: RNA, proteins, lipids, dyes, detergents, organic inhibitors of enzymatic reactions, buffers, salts, divalent**

Tissues and cells are desintegrated in the presence of lysis buffers combination. Further, Proteinase K digests cellular proteins, including DNA binding proteins, among them nucleases. Optimized buffer and ethanol are added to provide selective conditions for DNA binding to Magnetic Beads. Traces of contaminants

remaining in the solution are efficiently removed in three wash steps. High-quality cellular DNA is then eluted in low salt buffer, e.g.: Tris-HCl, TE or water. Isolated DNA is ready for downstream applications without the need for ethanol precipitation.

- **GeneMAGNET line is based on the use of silica paramagnetic beads (Magnetic Beads) for selective binding of RNA and DNA. The use of specially designed binding and washing buffers enables the efficient purification of highly pure nucleic acids.**



EURx Ltd. 80-297 Gdansk Poland  
ul. Przyrodnikow 3, NIP 957-07-05-191  
KRS 0000202039, [www.eurx.com.pl](http://www.eurx.com.pl)  
orders: email: [orders@eurx.com.pl](mailto:orders@eurx.com.pl)  
tel. +48 58 524 06 97, fax +48 58 341 74 23

