

## GeneMATRIX DNA/RNA Extracol Kit

Phenol-based reagent for the isolation of total RNA, DNA and protein from cells and tissues in a set with minicolumns.

● **Cat. no. E3750**

EURx Ltd. 80-297 Gdansk Poland  
ul. Przyrodnikow 3, NIP 957-07-05-191  
KRS 0000202039, [www.eurx.com.pl](http://www.eurx.com.pl)  
orders: email: [orders@eurx.com.pl](mailto:orders@eurx.com.pl)  
tel. +48 58 524 06 97, fax +48 58 341 74 23



# Table of Contents

<b>Introductory Notes.....</b>	<b>3</b>
<b>Equipment and reagents to be supplied by user .....</b>	<b>4</b>
<b>DNA/RNA isolation .....</b>	<b>4</b>
<b>Part I Disruption and sample lysis.....</b>	<b>4</b>
<b>Part II Phase Separation .....</b>	<b>5</b>
<b>Part III DNA/RNA Precipitation .....</b>	<b>6</b>
<b>Part IV DNA Purification.....</b>	<b>6</b>
<b>Part V RNA Purification.....</b>	<b>7</b>
<b>Part VI Protein precipitation.....</b>	<b>8</b>
<b>Safety Information .....</b>	<b>10</b>

<b>Content</b>	<b>25 preps E3750-01</b>	<b>100 preps E3750-02</b>	<b>Storage/Stability</b>
DNA/RNA Extracol*	30 ml	120 ml	2-25°C
Wash RBW	54 ml	216 ml	15-25°C
PLB	4.5 ml	18 ml	15-25°C
Bromophenol Blue 0.5%	0.1 ml	0.4 ml	15-25°C
RNase-free water	3.6 ml	15 ml	15-25°C
Elution	3.6 ml	15 ml	15-25°C
DNA Binding Columns	25	2 x 50	2-25°C
RNA Binding Columns	25	2 x 50	2-25°C
Protocol	1	1	

\* Contains phenol (toxic and corrosive) and guanidine isothiocyanate (irritant) and may be a health hazard if not handled properly. Always work with DNA/RNA Extracol in a fume hold and always wear a lab coat, gloves and safety glasses. It can be stored both at room temperature and at 2-8°C (fridge).

# Introductory Notes

**NOTE 1 • Application.** DNA/RNA Extracol is a reagent for the isolation of total RNA and DNA from cell and tissue samples of human, animal, plant, yeast or bacterial origin. **DNA/RNA Extracol** is a monophasic solution of phenol, chaotropic salts and other components design to facilitate the isolation of nucleic acids.

**NOTE 2 • Maximum Sample Volume.** 1.0 ml of **DNA/RNA Extracol** is sufficient to isolate DNA and RNA from a maximum of 100 mg tissue or is sufficient to lyse  $5-10 \times 10^6$  animal, plant or yeast cells or  $1 \times 10^7$  bacterial cells. The sample volume should not exceed 10% of the volume of **DNA/RNA Extracol** used for homogenization. When isolating nucleic acids from human leukocytes, starting volume of blood should not exceed 1.5 ml per 1 ml **DNA/RNA Extracol** used for leukocytes lysis. The maximum column binding capacity for DNA is 25  $\mu\text{g}$ . Loading more than 25  $\mu\text{g}$  DNA may lead to DNA contamination of the RNA eluate. The maximum column binding capacity for RNA is 125  $\mu\text{g}$ . The maximum volume of the column reservoir is 650  $\mu\text{l}$ .

**NOTE 3 • Homogenization and lysis.** Efficient disruption and homogenization of the starting material is requirement for most kind of samples. It can be carried out directly in **DNA/RNA Extracol** or in **RL** buffer (E0310). When using **RL** buffer, volume of homogenized tissue should not exceed 10% of the volume of **DNA/RNA Extracol** used for DNA/RNA isolation.

**NOTE 4 • Sample Storage.** After the cells or tissues have been homogenized or lysed in **DNA/RNA Extracol**, samples can be sotred at  $-80^\circ\text{C}$  for at least one month.

**NOTE 5 • Kit Compounds Storage.** Keep all solutions tightly closed to avoid evaporation, resulting in components concentration changes. **DNA/RNA Extracol** solution and other buffers are stable at temperature  $2-25^\circ\text{C}$ . For long term storage, store all components of the kit at  $2-8^\circ\text{C}$ .

**NOTE 6 • Caution.** **DNA/RNA Extracol** contains phenol (toxic and corrosive) and guanidine isothiocyanate (irritant) and may be a health hazard if not handled properly. Always work with **DNA/RNA Extracol** in a fume hold and always wear a lab coat, gloves and safety glasses.

**NOTE 7 • Elution buffer is a low salt solution, that contains no metal ion chealators (e.g. EDTA) that can inhibit subsequent enzymatic reactions.** Elution buffer composition is suitable for downstream applications such as digestion with restriction enzymes, phosphorylation, ligation, Sanger sequencing, NGS etc. It is also possible to elute the DNA with Tris-HCl, water or TE.

The method is based on the extraction of nucleic acids aqueous solutions using organic solvents. After homogenizing the sample with **DNA/RNA Extracol**, chloroform (or 1-bromo-3-chloropropane) is added, and the homogenate is allowed to separate into a clear upper aqueous layer, an interphase, and a lower organic layer. Separation of nucleic acids between the phases is pH dependent. At pH higher than 6 RNA and DNA remains in the aqueous phase. The highly effective RNase inhibitory property of **DNA/RNA Extracol** protects the integrity of the DNA/RNA during lysis and results in the isolation of high-quality material. Nucleic acids can be precipitated from the aqueous layer with isopropanol or separated into DNA and RNA fraction using minicolumns and wash buffers.

## Equipment and reagents to be supplied by user

- Chloroform or 1-bromo-3-chloropropane, isopropanol, ethanol 75%, ethanol 96-100%.
- Optional: RL buffer (E0310) for sample homogenization.
- Optional: Lyse RBC buffer (E0326) for erythrocytes lysis.
- Refrigerated laboratory centrifuge or microcentrifuge, disposable gloves, sterile RNase-free pipet tips, sterile RNase-free 1.5-2 ml tubes, vorteks, equipment for sample disruption and homogenization.
- For protein precipitation from tissue, cell culture, blood (leukocytes): GE solution – 0.3 M guanidine hydrochloride in 95% ethanol, either  $\beta$ -mercaptoethanol (14.3 M,  $\beta$ -ME) or [1 M] Dithiothreitol (DTT).
- For protein precipitation from plant tissue: either  $\beta$ -mercaptoethanol (14.3 M,  $\beta$ -ME) or [1 M] Dithiothreitol (DTT).

# DNA/RNA isolation

## Part I Disruption and sample lysis

### 1. Tissue:

Homogenize tissue samples in 1 ml of **DNA/RNA Extracol** per 10-100 mg of tissue. The tissue can also be homogenized in **RL** buffer (E0310). When using **RL** buffer for homogenization, in the next step, suspension of homogenized tissue must be added to **DNA/RNA Extracol** solution. Volume of homogenized tissue should not exceed 10% of the volume of **DNA/RNA Extracol** used for DNA/RNA isolation. For samples of fat tissue, a layer of fat may accumulate at the top, which should be removed. Centrifuge sample at 12 000 x g for 10 min at 4°C. Remove and discard the fatty layer.

### 2. Plant tissue:

Homogenize plant tissue samples in 1 ml of **DNA/RNA Extracol**. The tissue can also be homogenized in **RL** buffer (E0310). When using **RL** buffer for homogenization, in the next step, suspension of homogenized tissue must be added to **DNA/RNA Extracol** solution. Volume of homogenized tissue should not exceed 10% of the volume of **DNA/RNA Extracol** used for DNA/RNA isolation. For samples with high content polysaccharides or extracellular material, an additional centrifugation step is required to remove insoluble material from the sample. Centrifuge sample at 12 000 x g for 10 min at 4°C. Transfer the cleared supernatant to a new tube.

### 3. Cell grown in suspension:

Pellet cells by centrifugation and remove media. Lyse cells with 1 ml of DNA/RNA Extracol per  $5\text{-}10 \times 10^6$  cells and pass the lysate several times through a pipette tip.

### 4. Cell grown on monolayer:

Remove growth media. Lyse cells directly in a culture dish or flask by adding 1 ml of **DNA/RNA Extracol** per  $10 \text{ cm}^2$  growth area. Pipette the cell lysate several times to ensure sufficient cell disruption.

o *DNA/RNA Extracol is not compatible with plastic culture plates.*

### 5. Blood (leukocytes):

**DNA/RNA Extracol** can be used for DNA/RNA isolation from leukocytes. The maximum amount of human blood is 1.5 ml per 1 ml **DNA/RNA Extracol**. If the main purpose is the isolation of RNA, do not use frozen blood.

Add 4 volumes of buffer **Lyse RBC** (E0326) to a fresh blood. Mix by inverting the tube. Keep at  $4^\circ\text{C}$  for 10 min to lyse erythrocytes. Mix twice by inverting the tube. Centrifuge at  $400 \times g$  for 10 min at  $4^\circ\text{C}$ , and carefully decant the supernatant.

Add **DNA/RNA Extracol** to the leukocytes pellet. Mix thoroughly by pipetting for homogenization.

## Part II Phase Separation

1. Incubate samples for 5 min at room temperature.

2. Add 0.2 ml of chloroform (or 0,1 ml 1-bromo-3-chloropropane) per 1 ml of **DNA/RNA Extracol** used for homogenization.

3. Cover the sample tightly, shake vigorously for 15 sec.

4. Incubate samples for 2-5 min at room temperature.

5. Centrifuge sample at  $12\ 000 \times g$  for 15 min at  $4^\circ\text{C}$ .

o *Centrifugation separates the mixture into 3 phases: orange organic phase (containing protein), an interphase and a colorless upper aqueous phase (containing DNA and RNA). The upper aqueous phase is ~50% of the total volume.*

6. Remove the aqueous phase very carefully, without disturbing the interphase. For isolation of both, DNA and total RNA, continue with part III of the protocol **DNA/RNA Precipitation**. For separation of DNA from RNA continue with the steps described in part IV of the protocol **DNA Purification**.

Optionally, keep the residue (interphase and the bottom phenolic layer) to isolate the proteins. For protein purification, go to **Part VI Protein Precipitation**.

## Part III DNA/RNA Precipitation

1. Add 0.5 ml of 100% isopropanol to the aqueous phase, per 1 ml of **DNA/RNA Extracol** used for homogenization.
2. Incubate samples for 10 min at room temperature.
3. Centrifuge sample at 12 000 x g for 10 min at 4°C.
4. Remove the supernatant. Wash the pellet once with 1 ml 75% ethanol per 1 ml of **DNA/RNA Extracol** used in the initial homogenization. Vortex the samples briefly.
5. Centrifuge sample at 10 000 x g for 5 min at 4°C.
6. Remove the supernatant. Vacuum or air-dry the pellet and dissolve in **RNase-free water**.
  - o *Sample contains DNA, and RNA.*

## Part IV DNA Purification

1. The aqueous phase from the last step of part II of the protocol (Phase Separation) transfer to the **DNA binding spin-column** placed in a 2 ml receiver tube. Centrifuge at 12 000 x g for 1 min. Use the flow-through for RNA purification.
  - o *To receive RNA only, keep the flow-through and continue with part V of the protocol (RNA Purification).*
2. Add 600 µl of **Wash RBW** buffer to the DNA binding spin column and centrifuge at 12 000 x g for 1 min.
3. Remove the spin-column, pour off supernatant and place back into the receiver tube.
4. Add 300 µl of **Wash RBW** buffer to the spin column and centrifuge at 12 000 x g for 2 min.
  - o *Be careful not to contaminate the sample while removing spin-column from receiver tube. Check whether column is completely dry. If not, pour off supernatant and place back spin-column into the receiver tube and spin down for additional 1 min.*
5. Place spin-column into new receiver tube (1.5-2 ml) and add 50-100 µl of **Elution** buffer directly onto the membrane to elute bound DNA. Incubate spin-column/receiver tube assembly for 2 min at room temperature.
  - o *In order to improve the efficiency of the elution genomic DNA from membrane, Elution buffer can be heated to a temperature of 80°C.*
6. Centrifuge for 1 min at 12 000 x g. Remove spin-column, cap the receiver tube. Genomic DNA is ready for analysis/manipulations. It can be stored either at 2-8°C (preferred) or at -20°C (avoid multiple freezing and defrosting of DNA).

## Part V RNA Purification

1. To the flow-through from the step 2 of part IV of the protocol (DNA Purification) add the same volume of ethanol (96-100% [v/v]). Mix thoroughly. Do not centrifuge.
  - For example, if 400  $\mu$ l flow-through was collected in step 1 in Part IV of the protocol, add 400  $\mu$ l ethanol.
2. Apply up to 600  $\mu$ l of the mixture to the **RNA binding spin-column** and centrifuge at 12 000 x g for 1 min. Remove the spin-column, pour off supernatant and place back into the receiver tube.
3. Transfer the remaining mixture to the same **RNA binding spin-column** and centrifuge at 12 000 x g for 1 min. Remove the spin-column, pour off supernatant and place back into the receiver tube.
4. Add 600  $\mu$ l of **Wash RBW** buffer to the RNA binding spin column and centrifuge at 12 000 x g for 1 min.
5. Remove the spin-column, pour off supernatant and place back into the receiver tube.
6. Add 300  $\mu$ l of **Wash RBW** buffer to the spin column and centrifuge at 12 000 x g for 2 min.
  - Be careful not to contaminate the sample while removing spin-column from receiver tube. Check whether column is completely dry. If not, pour off supernatant and place back spin-column into the receiver tube and spin down for additional 1 min.
7. Place spin-column into new receiver tube (1.5-2 ml) and add 50-100  $\mu$ l of **RNase-free water** directly onto the membrane to elute bound RNA. Incubate spin-column/receiver tube assembly for 2 min at room temperature.
8. Centrifuge for 1 min at 12 000 x g. Remove spin-column, cap the receiver tube. RNA is ready for analysis/manipulations. Store the samples at -20°C or below.

## Part VI Protein precipitation

### Tissue, cell culture, blood (leukocytes)

1. To the phenol phase (step 6 in part II of the protocol), add 1.5 volumes of isopropanol. Mix thoroughly. Incubate for 10-30 min at room temperature.
2. Centrifuge at 12 000 x g for 10 min at 4°C, and carefully decant the supernatant.
3. To the protein pellet add 1 ml of GE solution, per 1 ml of DNA/RNA Extracol used for homogenization. Mix thoroughly. Incubate for 10-20 min at room temperature.
  - Prepare a wash solution GE consisting of 0.3 M guanidine hydrochloride in 95% ethyl alcohol.
4. Centrifuge at 7 500 x g for 5 min at 4°C, and carefully decant the supernatant.
5. Repeat steps 3 and 4.
6. Wash the pellet once with 2 ml of ethanol (96-100%). Vortex the samples briefly.
7. Centrifuge at 7 500 x g for 5 min at 4°C, and carefully decant the supernatant.
8. Dry the protein pellet for 5-15 min at room temperature.
9. Dissolve the protein pellet in 80-150 µl protein loading buffer **PLB** (Note 2).
  - Buffer PLB is a sample buffer for use in SDS-PAGE analysis. If the proteins will not be analyzed by SDS-PAGE, use a buffer compatible with the intended application. As a result of the method of isolation the precipitated protein is highly denatured and shows reduced solubility in water. Dissolution the precipitate is possible in PLB buffer or other solution containing a high concentration of detergent (eg 1.5-5% SDS). Therefore, Bradford and Lowry assays are not applicable for quantifying protein yield. For protein quantitation, use the Bicinchoninic Acid Assay (BCA).
  - For SDS-PAGE analysis, add 10 µl Bromophenol Blue and either 25 µl β-mercaptoethanol (β-ME) or 25 µl [1 M] dithiothreitol per 1 ml buffer PLB before use. After addition of β-ME store buffer PLB at 2-8°C. DTT-supplemented buffer PLB must not be stored, always prepare fresh aliquots.
  - In case of PLB buffer ingredients precipitation warm up until clarified.
10. Incubate for 5-10 min at 60°C to dissolve and denature sample.
11. If some insoluble material is still visible, centrifuge at maximum speed for 1 min. The supernatant is ready to use in downstream applications such as SDS-PAGE and others.
  - Sample can be stored at 2-8°C for short period or at -20°C for several months.



## Plant tissue

1. To the phenol phase (step 6 in part II of the protocol), add 2 volumes of ethanol (96-100%). Mix thoroughly. Incubate for 30 min at 4°C.
2. Centrifuge at 12 000 x g for 10 min at 4°C, and carefully decant the supernatant.
3. Wash the pellet with 2 ml of 75% ethanol per 1 ml of **DNA/RNA Extracol** used in the initial homogenization. Vortex the samples briefly.
4. Centrifuge at 7 500 x g for 5 min at 4°C, and carefully decant the supernatant.
5. Repeat steps 3 and 4.
6. Dry the protein pellet for 5-15 min at room temperature.
7. Dissolve the protein pellet in 80-150 µl protein loading buffer **PLB** (Note 2).
  - Buffer PLB is a sample buffer for use in SDS-PAGE analysis. If the proteins will not be analyzed by SDS-PAGE, use a buffer compatible with the intended application. As a result of the method of isolation the precipitated protein is highly denatured and shows reduced solubility in water. Dissolution the precipitate is possible in PLB buffer or other solution containing a high concentration of detergent (eg 1.5-5% SDS). Therefore, Bradford and Lowry assays are not applicable for quantifying protein yield. For protein quantitation, use the Bicinchoninic Acid Assay (BCA).
  - For SDS-PAGE analysis, add 10 µl Bromophenol Blue and either 25 µl β-mercaptoethanol (β-ME) or 25 µl [1 M] dithiothreitol per 1 ml buffer PLB before use. After addition of β-ME store buffer PLB at 2-8°C. DTT-supplemented buffer PLB must not be stored, always prepare fresh aliquots.
  - In case of PLB buffer ingredients precipitation warm up until clarified.
8. Incubate for 5-10 min at 60°C to dissolve and denature sample.
9. If some insoluble material is still visible, centrifuge at maximum speed for 1 min. The supernatant is ready to use in downstream applications such as SDS-PAGE and others.
  - Sample can be stored at 2-8°C for short period or at -20°C for several months.

# Safety Information

## DNA/RNA Extracol

### Danger



**H301+H311+H331** Toxic if swallowed, in contact with skin or if inhaled.

**H314** Causes severe skin burns and eye damage.



**H341** Suspected of causing genetic defects.  
**H373** May cause damage to organs through prolonged or repeated exposure.



**H411** Toxic to aquatic life with long lasting effects.



**P273** Avoid release to the environment.  
**P280** Wear protective gloves/protective clothing/eye protection/face protection.

**P301+P310** If swallowed: immediately call a poison center/doctor.

**P302+P352** If on skin: wash with plenty of water.

**P304+P340** If inhaled: remove person to fresh air and keep comfortable for breathing.

**P305+P351+P338** If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**P405** Store locked up.

**EUH032** Contact with acids liberates very toxic gas.

## Wash RBW

### Danger



**H225** Highly flammable liquid and vapour.

**H319** Causes serious eye irritation.

**P210** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.



**P280** Wear protective gloves/protective clothing/eye protection/face protection.

**P305+P351+P338** If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**P403+P235** Store in a well-ventilated place. Keep cool.

**P337+P313** If eye irritation persists: Get medical advice/ attention.



EURx Ltd. 80-297 Gdansk Poland  
ul. Przyrodnikow 3, NIP 957-07-05-191  
KRS 0000202039, [www.eurx.com.pl](http://www.eurx.com.pl)  
orders: email: [orders@eurx.com.pl](mailto:orders@eurx.com.pl)  
tel. +48 58 524 06 97, fax +48 58 341 74 23

