



Basic Cell Lysis

Cat. No.	size
E0294-01	50 ml
E0294-02	200 ml

Storage and transport: 4°C.

Term of validity: 12 months

Content: 10 mM Tris-HCl pH 7.4, 100 mM NaCl, 20 mM Na₄P₂O₇, 10 % glycerol, 1 % detergent, 2 mM NaF, 2 mM Na₃VO₄, 1 mM EDTA, 1 mM EGTA, 0.1 % SDS, 0.5 % deoxycholate

Product Important Information: Basic Cell Lysis does not contain protease inhibitors (1 mM PMSF or Protease Inhibitor Cocktail), which should be added directly before proteins extraction.

For extraction of nucleic and mitochondrial proteins it is recommended to use RIPA buffer (E0295).

Basic Cell Lysis contains phosphate inhibitors and is designed to lyse protein from adherent and suspension cultured mammalian cells and tissues, depending on nonionic detergent and recommended for cytoplasmic and membrane-bound proteins extraction.

The buffer is complemented by NaF and Na₃VO₄ and can be used to study the phosphorylated proteins.

Procedure for Lysis of Monolayer-cultured Adherent Mammalian Cells:

- Carefully remove culture medium from adherent cells.
- Wash cells twice with chilled PBS.
- Add cold Basic Cell Lysis to the cells (1 mL of buffer per 10⁸ cells - the volume might need optimization for different cell types and amount of the protein of interest). Incubate on ice for 15 minutes, swirling the plate occasionally until the clear lysate is visible.
- Collect the lysate and centrifuge at ~14,000 × g for 15 minutes.
- Transfer supernatant to a new tube for further analysis or store at -80 °C.

Procedure for Lysis of Suspension-cultured Mammalian Cells:

- Pellet the cells by centrifugation at 2500 × g for 5 minutes. Carefully remove and discard the supernatant.
- Wash cells twice in chilled PBS. Centrifuge at 2500 × g for 5 minutes. Carefully remove and discard the supernatant.
- Add cold Basic Cell Lysis to the cell pellet (1 mL of buffer per 10⁸ cells - the volume might need optimization for different cell types and amount of the protein of interest). Pipette the mixture up and down to suspend the pellet.
- Shake mixture gently for 15 minutes on ice until the clear lysate is visible.
- Centrifuge samples at 14,000 × g for 15 minutes.
- Transfer supernatant to a new tube for further analysis or store at -80 °C.

Procedure for Lysis of Tissues:

- Place the fresh tissue into chilled PBS and rinse several times. Mince the tissue into small pieces.
- Add cold Basic Cell Lysis to the tissue at 10:1 (10 mL Basic Cell Lysis per gram of tissue).
- Homogenize for several minutes at high speed until no tissue chunks remain.
- Incubate on ice for 3 minutes.
- Centrifuge at 10000 × g for 10 minutes
- Transfer supernatant to a new tube for further analysis or store at -80 °C.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23