



## RIPA Lysis Buffer

Cat. No.	size
E0295-01	50 ml
E0295-02	200 ml

**Storage condition:** Store at 4°C.

**Contents:** 25mM Tris-HCl pH 7.6, 150mM NaCl, 1% detergent, 1% sodium deoxycholate, 0.1% SDS.

### Important Product Information:

RIPA Buffer does not contain protease or phosphatase inhibitors. If necessary use protease and phosphatase inhibitors immediately before use.

Some protein kinases and other enzymes may be sensitive to the components of the RIPA Buffer, resulting in their decreased activity. In such cases, use milder buffers: **Basic Cell Lysis (E0294)**, **Mild Cell Lysis (E0292)** or **Mild Cell Lysis PLUS (E0293)**.

**RIPA Buffer (Radioimmunoprecipitation Assay) is a complete lysis solution reagent used for rapid and efficient total lysis and solubilization of protein from both adherent and suspension cultured mammalian cells, effectively extracting cytoplasmic, nuclear and membrane proteins.**

### Procedure for Lysis of Monolayer-cultured Adherent Mammalian Cells

- Carefully remove culture medium from adherent cells.
- Wash cells twice with chilled PBS.
- Add cold RIPA Buffer to the cells (1mL of buffer per  $5 \times 10^6$  HeLa or A431 cells). Incubate on ice for 5 minutes, swirling the plate occasionally.
- Collect the lysate and centrifuge at  $\sim 14,000 \times g$  for 15 minutes.
- Transfer supernatant to a new tube for further analysis.

### Procedure for Lysis of Suspension-cultured Mammalian Cells

- Pellet the cells by centrifugation at  $2500 \times g$  for 5 minutes. Carefully remove and discard the supernatant.
- Wash cells twice in chilled PBS. Centrifuge at  $2500 \times g$  for 5 minutes. Carefully remove and discard the supernatant.
- Add chilled RIPA Buffer to the cell pellet (1mL of RIPA buffer for  $5 \times 10^6$  of HeLa cells). Pipette the mixture up and down to suspend the pellet.
- Shake mixture gently for 15 minutes on ice. Centrifuge samples at  $14,000 \times g$  for 15 minutes.
- Transfer supernatant to a new tube for further analysis.

### Procedure for Lysis of Tissues

- Place the fresh tissue into chilled PBS and rinse several times. Mince the tissue into small pieces.
- Add RIPA Buffer to the tissue at 10:1 (10mL chilled RIPA Buffer per gram of tissue).
- Homogenize for several minutes at high speed until no tissue chunks remain.
- Incubate on ice for 3 minutes.
- Centrifuge at  $10000 \times g$  for 10 minutes
- Transfer supernatant to a new tube for further analysis.

### Troubleshooting

Problem	Possible Cause	Solution
Low total protein yield	Some cells are more resistant to lysis than others	Make sure the pellet is thoroughly suspended in RIPA Buffer and incubate for longer with occasional swirling – sonicate the pellet to increase yield
Low concentration of protein	Excess buffer used	Use less buffer
Proteolysis	No protease inhibitors added	Add protease inhibitor to the buffer before use
Low phosphorylation of proteins	Phosphatase activity	Add phosphatase inhibitor to the buffer before use

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039  
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23