



SG OneStep RT-qPCR kit

SG OneStep RT-qPCR kit is one-step RT-qPCR kit that provides accurate real-time quantification of RNA targets. Kit is composed of unique reverse transcriptase and highly processive hot start onTaq DNA Polymerase in easy to use format.

SG OneStep RT-qPCR kit

Kit Components

| Component | Cat. No. E0810-01 25 reactions of 25 µl | Cat. No. E0810-02 100 reactions of 25 µl |
|---|---|--|
| 2 x RT-qPCR SG Buffer | 1 x 350 µl | 2 x 0.7 ml |
| SG Enzyme Mix | 25 µl | 100 µl |
| Thermolabile UNG (uracil-N-glycosylase) 1 U/µl | 10 µl | 30 µl |
| Water, nuclease free | 1 x 0.5 ml | 2 x 1 ml |

SG OneStep RT-qPCR kit, plus ROX Solution

Kit Components

| Component | Cat. No. E0811-01 25 reactions of 25 µl | Cat. No. E0811-02 100 reactions of 25 µl |
|---|---|--|
| 2 x RT-qPCR SG Buffer | 1 x 350 µl | 2 x 0.7 ml |
| ROX Solution, 25 µM | 15 µl | 60 µl |
| SG Enzyme Mix | 25 µl | 100 µl |
| Thermolabile UNG (uracil-N-glycosylase) 1 U/µl | 10 µl | 30 µl |
| Water, nuclease free | 1 x 0.5 ml | 2 x 1 ml |

Storage

Store at -20°C in the dark.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23

The kit provides:

- 2 x RT-qPCR SG Buffer that is a universal reaction buffer with dNTPs (dTTP is partially replaced with dUTP) that can be used on most real-time PCR cyclers available. Allows as option to use thermolabile uracil-N-glycosylase (UNG).
- The SG Enzyme Mix contains unique highly sensitive reverse transcriptase, hot start onTaq DNA Polymerase, RNase Inhibitor and SYBR Green I dye.
- Reverse transcriptase that works in high range of temperatures from 35-55°C without loss of specificity and sensitivity.
- Both cDNA synthesis and PCR are performed in a single tube using gene-specific primers and either total RNA or mRNA.
- onTaq DNA Polymerase is a modified „hot start” enzymes which provides very tight inhibition of the polymerase activity at moderate temperatures which is restored during the initial denaturation step at 95°C for at least ten minutes.
- SYBR Green I is a fluorescent dye which binds all double-stranded DNA molecules and emits a fluorescent signal on binding. The excitation and emission maxima of SYBR Green I are at 494 nm and 521 nm, respectively, which are compatible with any real-time cycler.
- Kit contains thermolabile uracil-N-glycosylase (UNG) that is optimized for RT-qPCR reactions.
- If cyclers from Applied Biosystems are used ROX passive reference dye is necessary. SG OneStep RT-qPCR kit is provided in two variants: without ROX and with ROX Solution provided separately. The table below shows recommended amount of ROX (25 µM) required for a specific PCR cycler.

Recommended amounts of ROX for a specific real-time PCR cycler

| Instrument | Amount of ROX per 25 µl reaction | Final ROX concentration |
|---|--------------------------------------|-------------------------|
| Applied Biosystems: 7300, 7900HT, StepOne, StepOnePlus, ABI PRISM 7000 and 7700 | 0.5 µl | 500 nM |
| Applied Biosystems: 7500, ViiA 7, Stratagene: Mx3000P, Mx3005P, Mx4000 | 0.5 µl 10 x diluted (in water) | 50 nM |
| PCR machines from other manufacturers: Bio-Rad, Roche, Corbett, Eppendorf, Cepheid, etc. | Not required | - |

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23

Protocol

Preparation of PCR reaction:

| Component | Volume/reaction | Final concentration |
|---|----------------------------------|---------------------|
| 2 x RT-qPCR SG Buffer | 12.5 µl | 1 x |
| Forward Primer | Variable | 0.4 µM |
| Reverse Primer | Variable | 0.4 µM |
| Template RNA | Variable | ≤500 ng |
| SG Enzyme Mix | 1 µl | 1 µl/ reaction |
| Optional: ROX Solution, 25 µM | 0.5 µl or 0.5 µl 10 x diluted | 500 nM 50 nM |
| Optional: Thermolabile UNG (uracil-N-glycosylase) 1 U/µl | 0.25 µl | 0.25 U/reaction |
| Water, nuclease free | To 25 µl | - |
| Total volume | 25 µl | - |

Notes:

1. Minimize thaw-freeze cycles of 2 x RT-qPCR SG Buffer, keep SG Enzyme Mix and ROX solution on ice and minimize light exposure during handling to avoid loss of fluorescent signal intensity.
2. Thaw and gently vortex before use 2 x RT-qPCR SG Buffer.
3. A reaction volume of 25 µl should be used with most real-time cyclers. Other reaction volumes may be used if recommended for a specific instrument.
4. Optimal amplicon length in real-time RT-PCR using SYBR Green I is 70-150 bp.
5. To avoid amplification from genomic DNA design exon-exon primers.
6. Set up RT-PCR reactions on ice to minimize RNA template degradation.
7. The RNA template (≤500 ng/reaction) should be added to the individual PCR tubes or wells containing the whole reaction mix. Centrifuge briefly before placing into a cycler. Check if there are no bubbles left, if yes, spin again.
8. Place the samples in the cycler and start the program.
9. Reverse transcriptase works in a wide range of temperatures 35-55°C. The recommended starting temperature for reverse transcription is 50°C. For individual experiment temperature might be changed.
10. Standard concentration of MgCl₂ in real-time RT-PCR reaction is 3 mM (as provided with the 1 x RT-qPCR SG Buffer) in most cases this concentration will produce optimal results. However, if a higher MgCl₂ concentration is required, prepare a 25 mM MgCl₂ stock solution and add to reaction.
11. A final primer concentration of 0.3-0.5 µM is usually optimal, but can be individually optimized in the range of 0.1 µM to 1 µM. The recommended starting concentration is 0.4 µM. Raising primer concentration may increase PCR efficiency, but negatively affect RT-PCR specificity. Optimal primer concentration depends on the individual reaction and the real-time PCR cycler used.
12. Readjust the threshold value for analysis of every run.
13. If using Bio-Rad iCycler iQ or MyiQ instruments collect well factors at the beginning of each experiment. Use an external well factor plate according to the manufacturer's recommendations. Well factors are used to compensate for any excitation or pipetting variations.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23

Thermal Cycling Conditions:

| Step | Temperature | Time | Number of Cycles |
|--------------------------------------|-------------|------------|------------------|
| Reverse Transcription | 50°C | 20 min | 1 |
| Initial Denaturation | 95°C | 15 min | 1 |
| Denaturation | 94°C | 15 s | 40-45 |
| Annealing/Extension/Data acquisition | 60°C | 60 s | |
| Cooling | 4°C | Indefinite | 1 |

Notes:

1. During reverse transcription step of 50°C also thermolabile uracil-N-glycosylase might be used. Do not use UNG from *E.coli*, UNG will degrade all newly synthesized cDNA.
2. Thermolabile UNG is inactivated at 50°C during RT.
3. Melting curve analysis should be performed to verify the specificity and identity of PCR products. Melting curve analysis is an analysis step built into the software of real-time cyclers. Melting curve data between 65°C and 95°C should be acquired.
4. Data acquisition should be performed during the extension step. To suppress fluorescence readings caused by the generation of primer-dimers an additional data acquisition step can be added to the protocol. It is possible when T_m of primer-dimers differs from T_m of the specific product (T_m are generated during melting curve analysis). The temperature of the data acquisition step should be above T_m of primer-dimers but approximately 3°C below the T_m of the specific product.
5. Always check the RT-PCR product specificity by gel electrophoresis when designing a new assay. Melting temperatures of the specific product and primer-dimers may overlap.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23