



## T4 RNA Ligase 1

Cat. No.	size
E1059-01	1 000 units

**Concentration:** 10 000 units/ml

**Unit Definition:** One unit is defined as the amount of enzyme required to convert 1 nmole of 5'-[<sup>32</sup>P]rA<sub>16</sub> into a phosphatase-resistant form in 30 minutes at 37°C.

**Quality Control:** All preparations are assayed for contaminating endonucleases, exonucleases, nonspecific RNases, single- and double-stranded DNase activities. Greater than 90% as determined by SDS-PAGE.

### References:

1. Uhlenbeck, O.C., Gumpert, R.I., T4 RNA ligase, *The Enzymes*, 15B, 31-60, Academic Press Inc., 1982.
2. Middleton, T., et al., Synthesis and purification of oligoribonucleotides using T4 RNA ligase and reverse-phase chromatography, *Anal. Biochem.*, 144, 110-117, 1985.
3. Brennan, C.A., et al., Using T4 RNA ligase with DNA substrates, *Meth. Enzymol.*, 100, 38-52, 1983.
4. Tessier, D.C., et al., Ligation of single-stranded oligodeoxyribonucleotides by T4 RNA ligase, *Anal. Biochem.*, 158, 171-178, 1986.
5. Heckler, T.G., et al., T4 RNA ligase mediated preparation of novel "chemically misacylated" tRNA P<sup>he</sup>s, *Biochemistry*, 23, 1468-1473, 1984.
6. Edwards, J.B., et al., Oligodeoxyribonucleotide ligation to single-stranded cDNAs: a new tool for cloning 5'-ends of mRNAs and for constructing cDNA libraries by in vitro amplification, *Nucleic Acids Res.*, 19, 5227-5232, 1991.
7. Kaluz, S., et al., Enzymatically produced composite primers: an application of T4 RNA ligase-coupled primers to PCR, *BioTechniques*, 19, 182-186, 1995.

T4 RNA Ligase 1 catalyzes the ligation of a 5' phosphoryl-terminated nucleic acid donor to a 3' hydroxyl-terminated nucleic acid acceptor through the formation of a 3' → 5' phosphodiester bond with hydrolysis of ATP to AMP and PPi. Substrates include single-stranded RNA and DNA as well as dinucleoside pyrophosphates.

### Applications:

- RNA 3' -end labeling with Cytidine 3',5' - bis [α-32P] phosphate (1).
- Joining RNA to RNA (2).
- Synthesis of oligoribonucleotides and oligodeoxyribonucleotides (3, 4).
- Specific modifications of tRNAs (5).
- Oligodeoxyribonucleotide ligation to single-stranded cDNAs for 5'-RACE (Rapid Amplification of cDNA Ends) (6).
- Site-specific generation of composite primers for PCR (7).

**Source:** An *E. coli* strain containing a clone of the T4 RNA ligase gene.

**Format:** liquid

**Storage Buffer:** 20 mM HEPES (pH 7.4), 1 mM EDTA, 2 mM DTT, 50% glycerol.

**T4 RNA Ligase Reaction Buffer (1x):** 50 mM Tris-HCl (pH 7.8), 10 mM MgCl<sub>2</sub>, 1 mM DTT.

### RNA-RNA ligation exemplary reaction (20 μl):

Component:	Amount:
10 x T4 RNA Ligase Reaction Buffer	2 μl
RNase Inhibitor (E4210)	20 U
ATP	1 mM
DMSO (optional)	10%
ssRNA acceptor with free 3' -OH	20 pmol
RNA donor with free 5' -PO <sub>4</sub> and blocked 3' end	up to 200 pmol
T4 RNA Ligase	10 U
RNase-free Water	to 20 μl

Incubate at 25°C for 2 hours or at 16°C for 16 hours, proceed to RNA clean up for further RNA manipulation steps.