

## OMNI Nuclease


Cat. No.	size
E1120-01	20 000 units
E1120-02	100 000 units

### Unit Definition:


One unit is defined as the amount of enzyme that causes a  $\Delta A_{260}$  of 1.0 in 30 min, which corresponds to complete digestion of 37  $\mu$ g of Herring sperm DNA.

### Storage Conditions:

Store at -20°C. **Do not store at -80°C.**



**OMNI Nuclease activity test (1% agarose gel).** 37  $\mu$ g of Herring sperm DNA (1-4) and Torula Yeast RNA type IV (5-8) was incubated for 30 min at 37°C with increasing amount of OMNI Nuclease units.



### References:

1. Nestle M., Roberts W.K., (1969) *J. Biol. Chem.* 244, 5219-5225
2. Friedhoff P., Gimadutdinow O., (1994) *Protein Expr.Purif.* 5, 37-43
3. Meiss G., Friedhoff P., (1995) *Biochemistry* 34, 11979-11988

OMNI Nuclease catalyzes the removal of all forms of DNA and RNA (linear, circular, double and single stranded) (1). Nuclease completely degrades nucleic acids to 5'-monophosphate terminated oligonucleotides of 2-6 in length (2, 3).

### Applications:

- Viscosity reduction in bacterial, yeast and mammalian protein extracts.
- Sample preparation for protein 2D electrophoresis.
- Removal of nucleic acids contaminants from protein preparations.

### Assay Conditions:

50 mM Tris-HCl (pH 8.0), 1 mM MgCl<sub>2</sub>, 100  $\mu$ g/ml BSA, 37  $\mu$ g of Herring sperm DNA. Incubation is carried out at 37°C for 30 min in a reaction volume of 30  $\mu$ l.

### Storage Buffer:

20 mM Tris-HCl (pH 8.0), 20 mM NaCl, 2 mM MgCl<sub>2</sub>, 50% glycerol.

### Quality Control:

All preparations are assayed for contaminating protease activity.

### Instructions for use:

#### For removal of nucleic acids:

The substrate DNA or RNA dilute in assay buffer (50 mM Tris-HCl pH 8.0, 1 mM MgCl<sub>2</sub>, 100  $\mu$ g/ml BSA) to 2  $\mu$ g/ $\mu$ l. Incubate the substrate with different units of OMNI Nuclease at 37°C for 30 min. Analyze by gel electrophoresis. To remove trace DNA/RNA use 2 U of nuclease to 1  $\mu$ g of nucleic acids and extend incubation time to 24 h.

#### For viscosity reduction:

It is recommended to use 25 U/ml of lysate. The ratio of lysis buffer (ml) to the gram of cells should be 1 g : 1-10 ml

#### Example reaction:

1. 5g of *E. coli* (wet weight) resuspend in 25 ml of Tris-HCl pH 8.0, 1 mM EDTA (lysozyme, protease inhibitors if needed).
2. Add MgCl<sub>2</sub> the final concentration of 5 mM.
3. Add 2  $\mu$ l (500 U) of OMNI Nuclease.
4. Use disruption method of your choice (Ultrasonic disruption, French press, etc).
5. Obtain clear cell lysate supernatant by centrifugation ~12,000 rpm for 30 min at 4°C.

#### Notes:

1. OMNI Nuclease is very stable, can be diluted with 20 mM Tris-HCl (pH 8.0 at 25°C), 20 mM NaCl, 2 mM MgCl<sub>2</sub> and stored at 4°C for many days without loss of activity.
2. Displays a broad pH tolerance (from 6 to 10) and has a wide temperature spectrum (0-44°C).
3. OMNI Nuclease is not recommended for purification of proteins that must be nuclease-free. However it can be separated from the target molecule using chromatography for example anion and cation exchange, hydrophobic interaction, hydroxyapatite or size exclusion. OMNI Nuclease does not bind to Ni-NTA resin.
4. OMNI Nuclease is inhibited (app. 50% of activity) by:
  - monovalent cation concentrations (Na<sup>+</sup>, K<sup>+</sup>, etc.) > 50 mM,
  - phosphate concentration > 20 mM,
  - ammonium sulphate concentration > 25 mM.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.