



TEV Protease

Cat. No.	size
E4310-01	1 000 units
E4310-02	10 000 units

Unit Definition:

One unit is the amount of enzyme required to cleave >85% of 3 µg of control substrate (35 kDa fusion protein) in 1 hour at 30°C.

Storage Conditions:

Store at -20°C. No autolysis. Protein is stable at least for 9 months.

References:

1. Polayes D.A. *Et al.* (1998) *Methods in Molecular Medicine* Vol.13: 169-183.
2. Dougherty, W.G. *et al.* (1989) *Virology* 172, 302.
3. Dougherty, W.G., and Parks, T.D. (1989) *Virology* 172, 145.
4. Dougherty, W.G. *et al.* (1988) *EMBO* 7,1281.
5. Kapust, R.B., *at al.* (2002a) *Biochem.Biophys. Res. Commun.*294: 949-955.
6. Nallamsrtty, S. *at al.* (2004) *Protein Expr Purif.* 38(1): 108-15.
7. Mohanty, A.K. *at al.* (2003) *Protein Expr Purif.* 27: 109-114.

Cysteine protease from Tobacco Etch Virus for the removal of affinity tags from fusion proteins under target protein friendly conditions. The recombinant protein contains a N-terminal Poly-His tag and a C-terminal Polyarginine tag for easy removal after cleavage. Genetically engineered for resistance against autolysis, as well as for improved activity and performance

Description:

- TEV protease is a catalytic part of the Nuclear Inclusion protein „a” (NIa) from tobacco etch virus (TEV).
- TEV is a cysteine protease that specifically recognizes and cleaves a linear epitope with general sequence E-X-X-Y-X-Q-(G/S)(2-5) (where X is any amino acid) (2-5). The cleavage occurs between Q and G/S. The most common sequence is ENLYFQG or ENLYFQS.
- Resistant to many widely used serine and cysteine protease inhibitors like: PMSF, AEBSF, TLCK, E-64, “Complete” protease inhibitor cocktail (Roche).
- Robust enzyme active in the wide range of different buffers (with NaCl varied from 0 to 0.4 M and in pH from 4 to 9, enzyme tolerates MES, acetate, phosphate, glycerol and sorbitol).
- Active in a broad temperature range from 4 to 30°C (the enzyme is 3 times less active at 4°C than at 30°C)(6).
- Sensitive to some detergents (7).
- Extremely useful for removing affinity tags from fusion proteins in conditions friendly for target protein.

1 x Reaction Buffer:

20 mM Tris-HCl (pH 8.0), 10 mM NaCl, 1 mM MgCl₂, 10 mM β-mercaptoethanol.

Storage Buffer:

0.4 M NaCl, 50 mM Tris-HCl (pH 7.5), 2 mM EDTA, 1 mM DTT, 50% (v/v) glycerol.

Quality Control:

Protease is greater than 95% single-band pure without non-specific protease contamination.

Cleavage is performed in 1 x TEV buffer for one hour at 30°C. Determination of optimal cleavage temperature (between 4°C and 30°C) and time are performed accordingly (1).