



EUR_x[®]

cennik / price list
01/2022

Jesteśmy razem od 1998

24

Together since 1998

1. PCR & REAL TIME PCR	2-7
THERMOSTABLE DNA POLYMERASES	2-4
PCR KITS	4-5
REAL-TIME PCR KITS	6
qPCR-BASED SPECIALIZED KITS	7
dNTPs/rNTPs	7
PROTEINS & CHEMICALS FOR PCR	7
2. TRANSCRIPTION & REVERSE TRANSCRIPTION	8-9
TRANSCRIPTION KITS	8
REVERSE TRANSCRIPTASES	8
REVERSE TRANSCRIPTION CHEMICALS & KITS	8-9
3. LAMP REACTIONS	9
4. MOLECULAR WEIGHT MARKERS	9-10
PROTEIN MOLECULAR WEIGHT MARKERS	9
DNA MOLECULAR WEIGHT MARKERS	10
FLUORESCENT-LABELED DNA MOLECULAR WEIGHT MARKERS	10
5. MESOTHERMAL POLYMERASES	10-11
6. NUCLEASES	11
7. LIGASES & CLONING KITS	12
8. SPECIAL PURPOSE ENZYMES & REAGENTS	12-13
DNA MODIFICATION	12
DNA REPLICATION	12
COT DNA	12
DNA - PRECIPITATION & VISUALISATION	13
ENZYME - STABILIZNG AGENTS	13
AGARASES	13
9. WATER, BUFFERS & SOLUTIONS	13-14
10. XELEX KITS	14
11. PROTEASES, ENZYMES & REAGENTS FOR PROTEOMICS	14
12. TISSUE CULTURE MEDIA COMPONENTS & REAGENTS	15
13. AGAROSSES	15
14. RESTRICTION ENDONUCLEASES	15-17
15. RNA PURIFICATION SYSTEMS	17-18
16. DNA PURIFICATION SYSTEMS	18-19
17. EMULSION KITS	19
18. CHEMICALS, BUFFERS & ACCESSORIES FOR DNA/RNA PURIFICATION SYSTEMS	20
19. COLUMNS FOR DNA/RNA PURIFICATION SYSTEMS	21

NEW PRODUCTS



1. Restriction Endonuclease Mlu I, reduced star activity	16
2. GeneMAGNET Human and Animal Tissue DNA Purification Kit	19
3. GeneMAGNET Food DNA Purification Kit	19
4. Lyse Buffer PVP	20
5. RNA Storage	20

PROMOTIONAL PRICES



1. HLA Taq DNA Polymerase	3
2. phi29 DNA Polymerase	11
3. MDA phi29 Kit	11

PRODUCT NAME & DESCRIPTION

CAT.NUMBER

PACKAGE

PRICE

1. PCR & REAL-TIME PCR

THERMOSTABLE DNA POLYMERASES

All thermostable DNA Polymerases are also available as Kits, which include a DNA Polymerase and a mix of carefully tested and calibrated Ultrapure dNTPs (catalog numbers indicated with letters EK).

Taq DNA Polymerases

Taq DNA Polymerase

B

E2500-01	200 U	€ 26.50
E2500-04	500 U	€ 58.00
E2500-02	1 000 U	€ 101.00
E2500-03	5 000 U	€ 500.00

Ultrapure, recombinant Taq DNA polymerase for general Use (5.0 U/μl).

Taq DNA Polymerase

EK2500-01	200 U	€ 33.50
EK2500-04	500 U	€ 75.50
EK2500-02	1 000 U	€ 136.00
EK2500-03	5 000 U	€ 675.00

Taq DNA Polymerase Kit containing polymerase and dNTPs.

Color Taq DNA Polymerase

E2510-01	200 U	€ 26.50
E2510-04	500 U	€ 58.00
E2510-02	1 000 U	€ 101.00
E2510-03	5 000 U	€ 500.00

The enzyme is supplemented with two inert gel tracking dyes (red and yellow), allowing direct loading of PCR products to agarose gels (1.0 U/μl).

Color Taq DNA Polymerase

EK2510-01	200 U	€ 33.50
EK2510-04	500 U	€ 75.50
EK2510-02	1 000 U	€ 136.00
EK2510-03	5 000 U	€ 675.00

Color Taq DNA Polymerase Kit containing polymerase and dNTPs.

Native Taq DNA Polymerase

E2504-01	200 U	€ 31.00
E2504-04	500 U	€ 69.00
E2504-02	1 000 U	€ 122.00
E2504-03	5 000 U	€ 600.00

Native Taq DNA polymerase for general Use (5.0 U/μl).

Native Taq DNA Polymerase

EK2504-01	200 U	€ 38.00
EK2504-04	500 U	€ 86.50
EK2504-02	1 000 U	€ 157.00
EK2504-03	5 000 U	€ 775.00

Native Taq DNA Polymerase Kit containing polymerase and dNTPs.

OptiTaq DNA Polymerases

OptiTaq DNA Polymerase

E2600-01	200 U	€ 31.00
E2600-04	500 U	€ 69.00
E2600-02	1 000 U	€ 122.00
E2600-03	5 000 U	€ 600.00

Proofreading blend of DNA polymerases for any template and long PCR (5.0 U/μl).

OptiTaq DNA Polymerase

EK2600-01	200 U	€ 38.00
EK2600-04	500 U	€ 86.50
EK2600-02	1 000 U	€ 157.00
EK2600-03	5 000 U	€ 775.00

OptiTaq DNA Polymerase Kit containing polymerase and dNTPs.

Color OptiTaq DNA Polymerase

E2610-01	200 U	€ 31.00
E2610-04	500 U	€ 69.00
E2610-02	1 000 U	€ 122.00
E2610-03	5 000 U	€ 600.00

The enzyme is supplemented with two inert gel tracking dyes (red and yellow), allowing direct loading of PCR products to agarose gels (1.0 U/μl).

Color OptiTaq DNA Polymerase

EK2610-01	200 U	€ 38.00
EK2610-04	500 U	€ 86.50
EK2610-02	1 000 U	€ 157.00
EK2610-03	5 000 U	€ 775.00

Color OptiTaq DNA Polymerase Kit containing polymerase and dNTPs.

Hot Start Taq DNA Polymerases

Perpetual Taq DNA Polymerase

B

E2700-01	200 U	€ 68.00
E2700-04	500 U	€ 154.00
E2700-02	1 000 U	€ 261.00
E2700-03	5 000 U	€ 1239.00

Monoclonal antibody automatic "hot start" PCR system (2.5 U/μl).

Perpetual Taq DNA Polymerase

EK2700-01	200 U	€ 75.00
EK2700-04	500 U	€ 171.50
EK2700-02	1 000 U	€ 296.00
EK2700-03	5 000 U	€ 1414.00

Perpetual Taq DNA Polymerase Kit containing polymerase and dNTPs.

Color Perpetual Taq DNA Polymerase

E2710-01	200 U	€ 68.00
E2710-04	500 U	€ 154.00
E2710-02	1 000 U	€ 261.00
E2710-03	5 000 U	€ 1239.00

The enzyme is supplemented with two inert gel tracking dyes (red and yellow), allowing direct loading of PCR products to agarose gels (1.0 U/μl).

Vert!

N new product


B bestseller

U unique

P promotion


PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
Color Perpetual Taq DNA Polymerase <i>Color Perpetual Taq DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK2710-01	200 U	€ 75.00
	EK2710-04	500 U	€ 171.50
	EK2710-02	1 000 U	€ 296.00
	EK2710-03	5 000 U	€ 1414.00
tiTaq DNA Polymerase <i>New generation "hot start" Taq DNA Polymerase (2.5 U/μl).</i>	E2715-01	200 U	€ 31.00
	E2715-04	500 U	€ 69.00
	E2715-02	1 000 U	€ 122.00
	E2715-03	5 000 U	€ 600.00
tiTaq DNA Polymerase <i>tiTaq DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK2715-01	200 U	€ 38.00
	EK2715-04	500 U	€ 86.50
	EK2715-02	1 000 U	€ 157.00
	EK2715-03	5 000 U	€ 775.00
onTaq DNA Polymerase <i>The modified "hot start" enzyme which is blocked at moderate temperatures and allows room temperature reaction setup (2.5 U/μl).</i>	E2713-01	200 U	€ 35.00
	E2713-04	500 U	€ 79.00
	E2713-02	1 000 U	€ 140.00
	E2713-03	5 000 U	€ 690.00
onTaq DNA Polymerase <i>onTaq DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK2713-01	200 U	€ 43.00
	EK2713-04	500 U	€ 98.00
	EK2713-02	1 000 U	€ 180.00
	EK2713-03	5 000 U	€ 890.00
HLA Taq DNA Polymerase <i>HLA Taq DNA Polymerase is a specially formulated enzyme for HLA SSP typing Using commercially available kits.</i>	E2717-01	200 U	€ 31.00
	E2717-04	500 U	€ 69.00
	E2717-02	1 000 U	€ 122.00
	E2717-03	5 000 U	€ 600.00
Hot Start OptiTaq DNA Polymerases			
Perpetual OptiTaq DNA Polymerase <i>Monoclonal antibody automatic "hot start" proofreading PCR system (2.5 U/μl).</i>	E2720-01	200 U	€ 79.00
	E2720-04	500 U	€ 179.00
	E2720-02	1 000 U	€ 313.00
	E2720-03	5 000 U	€ 1488.00
Perpetual OptiTaq DNA Polymerase <i>Perpetual OptiTaq DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK2720-01	200 U	€ 86.00
	EK2720-04	500 U	€ 196.50
	EK2720-02	1 000 U	€ 348.00
tiOptiTaq DNA Polymerase <i>Proofreading blend of DNA polymerases for any template and long PCR, new generation "hot start" technology (2.5 U/μl).</i>	E2725-01	200 U	€ 40.00
	E2725-04	500 U	€ 90.00
	E2725-02	1 000 U	€ 160.00
tiOptiTaq DNA Polymerase <i>tiOptiTaq DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK2725-01	200 U	€ 47.00
	EK2725-04	500 U	€ 107.50
	EK2725-02	1 000 U	€ 195.00
	EK2725-03	5 000 U	€ 930.00
Amplus DNA Polymerase			
Amplus DNA Polymerase <i>Amplus DNA Polymerase Kit containing polymerase (5 U/μl) and dNTPs.</i>	EK2900-01	100 U	€ 19.00
	EK2900-02	500 U	€ 86.50
tiAmplus DNA Polymerase <i>tiAmplus DNA Polymerase Kit containing polymerase (2.5 U/μl) and dNTPs.</i>	EK2930-01	100 U	€ 26.00
	EK2930-02	500 U	€ 120.00
Hybrid DNA Polymerase			
Hybrid DNA Polymerase <i>Extremely thermostable DNA polymerase, enabling efficient high fidelity PCR of genomic targets up to 12 kb and episomal targets Up to 20 kb; (2.0 U/μl).</i>	E2950-01	100 U	€ 27.00
	E2950-02	500 U	€ 112.00
Hybrid DNA Polymerase <i>Hybrid DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK2950-01	100 U	€ 29.00
	EK2950-02	500 U	€ 123.00
onHybrid DNA Polymerase <i>Extremely thermostable DNA polymerase, enabling efficient high fidelity PCR of genomic targets up to 12 kb and episomal targets Up to 20 kb; new generation "hot start" technology (2.0 U/μl).</i>	E2960-01	100 U	€ 37.00
	E2960-02	500 U	€ 153.00

Vertel





PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
onHybrid DNA Polymerase <i>onHybrid DNA Polymerase Kit containing polymerase and dNTPs, new generation "hot start" technology.</i>	EK2960-01 EK2960-02	100 U 500 U	€ 44.00 € 185.00
Pfu DNA Polymerases			
Pfu DNA Polymerase <i>Highly thermostable DNA polymerase for PCR applications requiring high fidelity, proofreading activity; (5.0 U/μl).</i>	E1114-01 E1114-02 E1114-03	100 U 500 U 2 500 U	€ 20.00 € 75.00 € 370.00
Pfu DNA Polymerase <i>Pfu DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK1114-01 EK1114-02 EK1114-03	100 U 500 U 2 500 U	€ 22.00 € 87.00 € 410.00
PfuPlus! DNA Polymerase <i>Extremely thermostable proofreading DNA polymerase blend, formulated for efficient site-directed mutagenesis and synthesis of DNA products Up to 20 kb (5.0 U/μl).</i>	E1118-01 E1118-02 E1118-03	100 U 500 U 2 500 U	€ 20.00 € 75.00 € 370.00
PfuPlus! DNA Polymerase <i>PfuPlus! DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK1118-01 EK1118-02 EK1118-03	100 U 500 U 2 500 U	€ 22.00 € 87.00 € 410.00
Color PfuPlus! DNA Polymerase <i>The enzyme is supplemented with two inert gel tracking dyes (red and yellow), allowing direct loading of PCR products to agarose gels (1.25 U/μl).</i>	E1110-01 E1110-02 E1110-03	100 U 500 U 2 500 U	€ 20.00 € 75.00 € 370.00
Color PfuPlus! DNA Polymerase <i>Color PfuPlus! DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK1110-01 EK1110-02 EK1110-03	100 U 500 U 2 500 U	€ 22.00 € 87.00 € 410.00
onPfuPlus! DNA Polymerase <i>Extremely thermostable proofreading DNA polymerase blend, formulated for efficient site-directed mutagenesis and synthesis of DNA products Up to 20 kb; new generation "hot start" technology (5.0 U/μl).</i>	E1113-01 E1113-02 E1113-03	100 U 500 U 2 500 U	€ 29.00 € 123.00 € 535.00
onPfuPlus! DNA Polymerase <i>onPfuPlus! DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK1113-01 EK1113-02 EK1113-03	100 U 500 U 2 500 U	€ 31.00 € 135.00 € 575.00
Thermus DNA Polymerases			
Tfi DNA Polymerase <i>Native Thermus flavus DNA polymerase (5.0 U/μl).</i>	E1112-01 E1112-02	200 U 1 000 U	€ 50.00 € 198.00
Tfi DNA Polymerase <i>Tfi DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK1112-01 EK1112-02	200 U 1 000 U	€ 57.00 € 233.00
Tth DNA Polymerase <i>Thermus thermophilus DNA polymerase (5.0 U/μl).</i>	E1115-01 E1115-02	100 U 500 U	€ 45.00 € 180.00
Tth DNA Polymerase <i>Tth DNA Polymerase Kit containing polymerase and dNTPs.</i>	EK1115-01 EK1115-02	100 U 500 U	€ 48.50 € 197.50
PCR KITS			
NXT Taq PCR Kit <i>Ready-to-use Solution containing "hot start" NXT Taq DNA Polymerase, reaction buffer, MgCl₂ and dNTPs. NXT Taq PCR Kit is designed for fast-cycling PCR on any thermal cycler.</i>	E2530-01 E2530-02 E2530-03	100 reactions 200 reactions 500 reactions	€ 40.00 € 78.00 € 160.00
NXT tiTaq PCR Kit <i>Ready-to-use Solution containing "hot start" NXT tiTaq DNA Polymerase, reaction buffer, MgCl₂ and dNTPs. NXT Taq PCR Kit is designed for fast-cycling PCR on any thermal cycler.</i>	E2535-01 E2535-02 E2535-03	100 reactions 200 reactions 500 reactions	€ 30.00 € 58.00 € 120.00
Taq PCR Master Mix (2x) <i>Ready-to-use Solution containing Taq DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2520-01 E2520-02 E2520-03	100 reactions 200 reactions 500 reactions	€ 32.00 € 60.00 € 130.00



Verte!



PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
Perpetual Taq PCR Master Mix (2x) <i>Ready-to-use Solution containing Perpetual Taq DNA Polymerase, antibodies, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2740-01	100 reactions	€ 75.00
	E2740-02	200 reactions	€ 145.00
	E2740-03	500 reactions	€ 290.00
Perpetual OptiTaq PCR Master Mix (2x) <i>Ready-to-use Solution containing Perpetual OptiTaq DNA Polymerase, antibodies, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2770-01	100 reactions	€ 75.00
	E2770-02	200 reactions	€ 145.00
	E2770-03	500 reactions	€ 290.00
OptiTaq PCR Master Mix (2x) <i>Ready-to-use Solution containing OptiTaq DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2910-01	100 reactions	€ 35.00
	E2910-02	200 reactions	€ 65.00
	E2910-03	500 reactions	€ 130.00
Color Taq PCR Master Mix (2x) <i>Ready-to-use Solution containing Taq DNA Polymerase, optimized and colored reaction buffer, MgCl₂ and dNTPs.</i>	E2525-01	100 reactions	€ 32.00
	E2525-02	200 reactions	€ 60.00
	E2525-03	500 reactions	€ 130.00
onTaq PCR Master Mix (2x) <i>Ready-to-use Solution containing onTaq DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2714-01	100 reactions	€ 45.00
	E2714-02	200 reactions	€ 87.00
	E2714-03	500 reactions	€ 180.00
tiTaq PCR Master Mix (2x) <i>Ready-to-use Solution containing tiTaq DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2716-01	100 reactions	€ 40.00
	E2716-02	200 reactions	€ 78.00
	E2716-03	500 reactions	€ 160.00
tiOptiTaq PCR Master Mix (2x) <i>Ready-to-use Solution containing tiOptiTaq DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2726-01	100 reactions	€ 45.00
	E2726-02	200 reactions	€ 87.00
	E2726-03	500 reactions	€ 180.00
Color Perpetual Taq PCR Master Mix (2x) <i>Ready-to-use Solution containing Perpetual Taq DNA Polymerase, antibodies, optimized and colored reaction buffer, MgCl₂ and dNTPs.</i>	E2745-01	100 reactions	€ 75.00
	E2745-02	200 reactions	€ 145.00
	E2745-03	500 reactions	€ 290.00
Hybrid PCR Master Mix (2x) <i>Ready-to-use Solution containing Hybrid DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2750-01	100 reactions	€ 40.00
	E2750-02	500 reactions	€ 160.00
onHybrid PCR Master Mix (2x) <i>Ready-to-use Solution containing onHybrid DNA Polymerase, optimized reaction buffer, MgCl₂ and dNTPs.</i>	E2755-01	100 reactions	€ 50.00
	E2755-02	500 reactions	€ 200.00
Color OptiTaq PCR Master Mix (2x) <i>Ready-to-use Solution containing OptiTaq DNA Polymerase, optimized and colored reaction buffer, MgCl₂ and dNTPs.</i>	E2915-01	100 reactions	€ 35.00
	E2915-02	200 reactions	€ 65.00
	E2915-03	500 reactions	€ 130.00
Multiplex PCR Master Mix (2x) <i>Ready-to-use Solution for multiplex PCR without the need for optimization.</i>	E2820-01	50 reactions	€ 60.00
	E2820-02	100 reactions	€ 110.00
Direct Blood PCR Kit <i>Enables to perform PCR reactions directly from whole blood with no prior DNA extraction or purification.</i>	E0950-01	100 reactions	€ 58.00
	E0950-02	500 reactions	€ 232.00
Direct Tissue PCR Kit <i>Direct Tissue PCR Kit enables to perform PCR reactions directly from tissue samples with no prior DNA purification.</i>	E0940-01	100 reactions	€ 46.00
	E0940-02	500 reactions	€ 210.00
Direct Plant PCR Kit <i>Enables to perform PCR reactions directly from plants with no prior DNA extraction or purification.</i>	E0960-01	100 reactions	€ 46.00
	E0960-02	500 reactions	€ 210.00

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
REAL-TIME PCR KITS			
SG qPCR Master Mix (2x) <i>For quantitative real-time PCR and two-step real-time RT-PCR.</i>	E0401-01 E0401-02 E0401-03	100 reactions 200 reactions 1 000 reactions	€ 43.00 € 82.00 € 370.00
SG qPCR Master Mix (2x), plus ROX Solution <i>For quantitative real-time PCR and two-step real-time RT-PCR. Contains ROX dye.</i>	E0402-01 E0402-02 E0402-03	100 reactions 200 reactions 1 000 reactions	€ 46.00 € 89.00 € 390.00
SG/ROX qPCR Master Mix (2x) <i>For quantitative real-time PCR and two-step real-time RT-PCR. Mix contains ROX dye.</i>	E0400-01 E0400-02 E0400-03	100 reactions 200 reactions 1 000 reactions	€ 46.00 € 89.00 € 390.00
Blank qPCR Master Mix (2x) <i>For quantitative real-time PCR and two-step real-time RT-PCR. The master mix does not contain any dye and allows the User to choose a dye and optimize the dye concentration.</i>	B E0403-01 E0403-02 E0403-03	100 reactions 200 reactions 1 000 reactions	€ 42.00 € 80.00 € 360.00
Blank qPCR Master Mix (2x), plus ROX Solution <i>For quantitative real-time PCR and two-step real-time RT-PCR. The master mix does not contain any dye and allows the User to choose a dye and optimize the dye concentration. Contains ROX dye.</i>	E0404-01 E0404-02 E0404-03	100 reactions 200 reactions 1 000 reactions	€ 45.00 € 87.00 € 380.00
SG onTaq qPCR Master Mix (2x) <i>For quantitative real-time PCR and two-step real-time RT-PCR, with onTaq DNA Polymerase.</i>	E0405-01 E0405-02 E0405-03	100 reactions 200 reactions 1 000 reactions	€ 42.00 € 80.00 € 360.00
SG onTaq qPCR Master Mix (2x), plus ROX Solution. <i>For quantitative real-time PCR and two-step real-time RT-PCR, with onTaq DNA Polymerase. Contains ROX dye.</i>	E0406-01 E0406-02 E0406-03	100 reactions 200 reactions 1 000 reactions	€ 45.00 € 87.00 € 380.00
SG/ROX onTaq qPCR Master Mix (2x) <i>For quantitative real-time PCR and two-step real-time RT-PCR. Mix contains onTaq polymerase and ROX dye.</i>	E0407-01 E0407-02 E0407-03	100 reactions 200 reactions 1 000 reactions	€ 45.00 € 87.00 € 380.00
Fast SG qPCR Master Mix (2x) <i>For fast-cycling quantitative real-time PCR and two-step real-time RT-PCR.</i>	E0411-01 E0411-02 E0411-03	100 reactions 200 reactions 1 000 reactions	€ 48.00 € 90.00 € 425.00
Fast SG qPCR Master Mix (2x), plus ROX Solution <i>For fast-cycling quantitative real-time PCR and two-step real-time RT-PCR. Contains ROX dye.</i>	E0412-01 E0412-02 E0412-03	100 reactions 200 reactions 1 000 reactions	€ 51.00 € 97.00 € 445.00
Fast SG/ROX qPCR Master Mix (2x) <i>For fast-cycling quantitative real-time PCR and two-step real-time RT-PCR. Mix contains ROX dye.</i>	E0413-01 E0413-02 E0413-03	100 reactions 200 reactions 1 000 reactions	€ 51.00 € 97.00 € 445.00
Probe qPCR Master Mix (2x) <i>For quantitative real-time PCR and two-step real-time RT-PCR Using sequence-specific probes.</i>	E0420-01 E0420-02 E0420-03	100 reactions 200 reactions 1 000 reactions	€ 46.00 € 89.00 € 390.00
Probe qPCR Master Mix (2x), plus ROX Solution <i>For quantitative real-time PCR and two-step real-time RT-PCR Using sequence-specific probes. Contains ROX dye.</i>	E0421-01 E0421-02 E0421-03	100 reactions 200 reactions 1 000 reactions	€ 49.00 € 96.00 € 410.00
Fast Probe qPCR Master Mix (2x) <i>For fast-cycling quantitative real-time PCR and two-step real-time RT-PCR Using sequence-specific probes.</i>	E0422-01 E0422-02 E0422-03	100 reactions 200 reactions 1 000 reactions	€ 46.00 € 89.00 € 390.00
Fast Probe qPCR Master Mix (2x), plus ROX Solution <i>For fast-cycling quantitative real-time PCR and two-step real-time RT-PCR Using sequence-specific probes. Contains ROX dye.</i>	E0423-01 E0423-02 E0423-03	100 reactions 200 reactions 1 000 reactions	€ 49.00 € 96.00 € 410.00
MP qPCR Master Mix (2x) <i>For quantitative multiplex real-time PCR and two-step real-time RT-PCR Using sequence-specific probes.</i>	E0424-01 E0424-02 E0424-03	100 reactions 200 reactions 1 000 reactions	€ 46.00 € 89.00 € 390.00
MP qPCR Master Mix (2x), plus ROX Solution <i>For quantitative multiplex real-time PCR and two-step real-time RT-PCR Using sequence-specific probes. Contains ROX dye.</i>	E0425-01 E0425-02 E0425-03	100 reactions 200 reactions 1 000 reactions	€ 49.00 € 96.00 € 410.00
HRM PCR Master Mix <i>For detection of gene mutations and SNPs by HRM analysis.</i>	E0440-01 E0440-02	100 reactions 500 reactions	€ 46.00 € 210.00



PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
qPCR-BASED SPECIALIZED KITS			
SARS-CoV-2 qPCR Detection Kit  For SARS -CoV-2 detection.	E0430-01 E0430-02	100 reactions 500 reactions	€ 540,00 € 2 750,00
N/RdRp SARS-CoV-2 Detection KIT For SARS-CoV-2 specific genes detection in samples of various origin	E0435-01 E0435-02	100 reactions 500 reactions	€ 500,00 € 2 450,00
SARS-CoV-2 Blue-LAMP Detection Kit  For SARS -CoV-2 detection (colorimetric analysis).	E0460-01 E0460-02	50 reactions 500 reactions	€ 200,00 € 1 800,00
Borrelia qPCR Detection Kit For Borrelia burgdorferi detection.	E0450-01 E0450-02	50 reactions 100 reactions	€ 150,00 € 290,00
Borrelia qPCR Detection Kit, plus ROX Solution For Borrelia burgdorferi detection. Contains ROX dye.	E0451-01 E0451-02	50 reactions 100 reactions	€ 152,00 € 295,00
dNTPs/rNTPs			
10 mM Ultrapure dNTPs Set Ultrapure set of dNTPs contains 4 separate 10 mM Solutions of dATP, dCTP, dGTP and dTTP.	E0501-01 E0501-02	4 x 200 µl 4 x 1 ml	€ 8,00 € 40,00
100 mM Ultrapure dNTPs Set Ultrapure set of dNTPs contains 4 separate 100 mM Solutions of dATP, dCTP, dGTP and dTTP.	E0502-01 E0502-02	4 x 200 µl 4 x 1 ml	€ 50,00 € 235,00
10 mM Ultrapure dNTPs Mix Mix of dNTPs contains dATP, dCTP, dGTP and dTTP in concentration of 10 mM each.	E0503-01 E0503-02	1 ml 5 x 1 ml	€ 50,00 € 200,00
5 mM Ultrapure dNTPs Mix Mix of dNTPs contains dATP, dCTP, dGTP and dTTP in concentration of 5 mM each.	E0505-01 E0505-02	200 µl 1 ml	€ 15,00 € 50,00
100 mM dUTP 100 mM Sodium salt Solution	E0515-01 E0515-02	200 µl 1 ml	€ 15,00 € 50,00
100 mM Ribo-Nucleotides Set Ultrapure set of rNTPs contains 4 separate 100 mM Solutions of rATP, rCTP, rGTP and rUTP.	E0601-01 E0601-02	4 x 200 µl 4 x 1 ml	€ 35,00 € 150,00
PROTEINS & CHEMICALS FOR PCR			
Gp 32 Protein Single-stranded DNA Binding Protein of T4 bacteriophage.	E1170-01 E1170-02	100 µg 500 µg	€ 54,00 € 225,00
Single-stranded DNA Binding Protein Single-stranded DNA Binding Protein from Escherichia coli.	E4200-01 E4200-02	100 µg 500 µg	€ 45,00 € 180,00
Taq Single-stranded DNA Binding Protein Single-stranded DNA Binding Protein from Thermus aquaticus.	E4300-01 E4300-02	50 µg 250 µg	€ 45,00 € 180,00
Uracil - N - Glycosylase (UNG) Eliminates Uracil from DNA molecules, prevents carryover contamination between reactions.	E1250-01	500 U	€ 49,00
Thermolabile Uracil - N - Glycosylase (UNG) Eliminates Uracil from DNA molecules. Heat inactivation (5 minutes at 50°C).	E1251-01	500 U	€ 49,00

PRODUCT NAME & DESCRIPTION

CAT.NUMBER

PACKAGE

PRICE

2. TRANSCRIPTION & REVERSE TRANSCRIPTION

TRANSCRIPTION KITS

T7 transcription Kit

T7 transcription Kit with modified T7 RNA Polymerase for higher tolerance towards modified nucleotides.

E0901-01

25 reactions

€ 100.00

E0901-02

50 reactions

€ 180.00

Apt-Get 2'-F T7 Transcription Kit 2'-Fluoro-Pyrimidine Set

T7 transcription Kit for synthesis of nuclease-resistant RNA

E0905-01

25 reactions

€ 100.00

E0905-02

50 reactions

€ 180.00

REVERSE TRANSCRIPTASES

AMV Reverse Transcriptase Native

Native RNA dependent DNA polymerase from Avian Myeloblastosis Virus.



E1372-01

500 U

€ 100.00

E1372-02

2 500 U

€ 475.00

HIV Reverse Transcriptase

RNA dependent DNA polymerase from Human Immunoeficiency Virus.

E1373-01

500 U

€ 270.00

E1373-02

2 500 U

€ 1080.00

Tth Reverse Transcriptase

RNA dependent DNA polymerase from Thermus thermophilus, active at elevated temperatures.

E1374-01

100 U

€ 72.00

E1374-02

500 U

€ 288.00

smART Reverse Transcriptase

Recombinant reverse transcriptase with reduced RNase H activity, increased thermalactivity and processivity.

E1376-01

10 000 U

€ 60.00

E1376-02

4 x 10 000 U

€ 215.00

REVERSE TRANSCRIPTION CHEMICALS & KITS

NG dART RT Kit

New Generation! For amplification of DNA from any RNA with high specificity and sensitivity in a two-step process.

E0801-01

25 reactions

€ 50.00

E0801-02

100 reactions

€ 165.00

E0801-03

50 reactions

€ 90.00

NG dART RT-PCR Kit

New Generation! For amplification of DNA from any RNA template with high specificity and sensitivity in a two-step process cDNA template is Used to amplify ds DNA of interest Using high fidelity OptiTaQ DNA Polymerase.

E0802-01

25 reactions

€ 55.00

E0802-02

100 reactions

€ 185.00

OneStep RT-PCR Kit

For amplification of DNA from any RNA with high specificity and sensitivity in single tube.

E0803-01

25 reactions

€ 60.00

E0803-02

100 reactions

€ 230.00

smART First Strand cDNA Synthesis Kit

cDNA synthesis Kit based on thermostabile and highly processive smART reverse transcriptase.

E0804-01

2 000 U

€ 30.00

E0804-02

10 000 U

€ 105.00

E0804-03

4 x 10 000 U

€ 340.00

smART RT-PCR Kit

Complete two-step RT-PCR Kit based on smART reverse transcriptase and Opti-TaQ DNA polymerase.

E0805-01

25 reactions

€ 60.00

E0805-02

100 reactions

€ 230.00

smART RT-qPCR Kit

Complete two-step RT-PCR Kit based on smART reverse transcriptase and SG (2x) qPCR master mix.

E0806-01

25 reactions

€ 60.00

E0806-02

100 reactions

€ 230.00

smART RT-qPCR Kit, plus ROX Solution

Complete two-step RT-PCR Kit based on smART reverse transcriptase and SG (2x) qPCR master mix with reference ROX dye.

E0807-01

25 reactions

€ 63.00

E0807-02

100 reactions

€ 236.00

Positive Control Set for RT-PCR

Set of primers and synthetic RNA template for positive control of RT-PCR reaction.

E0820-01

20 reactions

€ 25.00

SG OneStep qRT-PCR Kit

Universal solution for quantitative real-time RT-PCR in single tube.

E0810-01

25 reactions

€ 25.00

E0810-02

100 reactions

€ 100.00

SG OneStep qRT-PCR Kit, plus ROX Solution

Universal solution for quantitative real-time RT-PCR in single tube with reference ROX dye.

E0811-01

25 reactions


€ 28.00

E0811-02

100 reactions

€ 106.00

Verte!



new product



bestseller



unique



promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
Probe OneStep RT-qPCR Kit <i>Universal solution for quantitative real-time RT-PCR with probes in single tube.</i>	E0812-01 E0812-02	25 reactions 100 reactions	€ 25.00 € 100.00
Probe OneStep RT-qPCR Kit, plus ROX Solution <i>Universal solution for quantitative real-time RT-PCR with probes in single tube with reference ROX dye.</i>	E0813-01 E0813-02	25 reactions 100 reactions	€ 28.00 € 106.00
Probe Multiplex OneStep RT-qPCR Kit <i>Provides accurate real-time quantification of multiple RNA targets in gene expression analysis with probes in single tube.</i>	E0819-01 E0819-02	100 reactions 500 reactions	€ 135.00 € 555.00
OneStep PRO RT-qPCR Kit <i>Fully thermostable and high fidelity system for quantitative real-time RT-PCR in single tube.</i>	E0825-01 E0825-02	25 reactions 100 reactions	€ 25.00 € 100.00
OneStep PRO RT-qPCR Kit, plus ROX Solution <i>Fully thermostable and high fidelity system for quantitative real-time RT-PCR in single tube with reference ROX dye.</i>	E0826-01 E0826-02	25 reactions 100 reactions	€ 28.00 € 127.00
Ribonuclease Inhibitor <i>Ribonuclease inhibitor suitable for Use in enzymatic reactions (RNase-free).</i>	B E4210-01 E4210-02	7 500 U 37 500 U	€ 95.00 € 380.00
Placental Ribonuclease Inhibitor <i>Recombinant human placental protein which specifically inhibits ribonucleases A, B and C.</i>	E4215-01 E4215-02	2 000 U 10 000 U	€ 48.00 € 190.00
Random Hexamer Primers <i>A mixture of single-stranded random hexonucleotides with 5' and 3' - hydroxyl ends (5' NNNNNN 3').</i>	E0101-01	20 µg	€ 25.00
Modified Random Hexamer Primers <i>A mixture of single-stranded random hexonucleotides with 5' and 3' - hydroxyl ends (5' NNNN*N*N 3'), resistant to exonuclease 3'->5' degradation. *-PTO</i>	E0107-01	100 µl	€ 50.00
Oligo(dT)₂₀ Primers <i>The solution of oligonucleotides composed of 20 deoxythymidylic acid residues.</i>	E0105-01	25 µg	€ 50.00
Anchored Oligo(dT)₂₀ VN Primers <i>The primer mixture of a string of 20 thymidine residues followed by dV at 3' end.</i>	E0106-01	25 µg	€ 65.00


3. LAMP REACTIONS

Blue-LAMP Kit <i>Isothermal DNA amplification kit with colorimetric detection</i>	E1410-01 E1410-02	100 reactions 500 reactions	€ 100.00 € 440.00
Blue-LAMP RT Kit <i>Isothermal DNA and RNA amplification kit with colorimetric detection</i>	E1415-01 E1415-02	100 reactions 500 reactions	€ 110.00 € 495.00
Flo-LAMP Kit <i>Isothermal DNA amplification kit with fluorescence detection</i>	E1420-01 E1420-02	100 reactions 500 reactions	€ 110.00 € 495.00
Flo-LAMP RT Kit <i>Isothermal DNA and RNA amplification kit with fluorescence detection</i>	E1425-01 E1425-02	100 reactions 500 reactions	€ 120.00 € 540.00

4. MOLECULAR WEIGHT MARKERS

PROTEIN MOLECULAR WEIGHT MARKERS

Perfect™ Color Protein Ladder <i>Prestained protein ladder for sizing polypeptides ranging from 7 kDa to 240 kDa (7, 15, 20, 25, 35, 50, 70, 100, 140, 240 kDa).</i>	E3215-01	50 loadings	€ 55.00
Perfect™ Tricolor Protein Ladder <i>Prestained protein ladder for sizing polypeptides ranging from 11 kDa to 245 kDa (11, 17, 20, 25, 35, 48, 63, 75, 100, 135, 180, 245 kDa).</i>	E3210-01	50 loadings	€ 55.00



new product



bestseller




unique



promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
DNA MOLECULAR WEIGHT MARKERS			
Perfect Plus™ 50-500 bp DNA Ladder <i>Ready to Use ladder for sizing small DNA fragments from 50 to 500 bp (50, 100, 250, 400, 500 bp).</i>	E3145-01 E3145-02	100 loadings 500 loadings	€ 55.00 € 220.00
Perfect Plus™ 2 kb DNA Ladder <i>Ready to Use DNA ladder for sizing DNA fragments from 50 - 2000 bp (50, 150, 300, 500, 750, 1000, 1500, 2000 bp).</i>	E3140-01 E3140-02	100 loadings 500 loadings	€ 55.00 € 220.00
Perfect Plus™ Molecular Weight Quantitative DNA <i>Ready to Use DNA ladder for sizing and quantifying DNA fragments from 100 to 5000 bp (100, 250, 500, 750, 1000, 1500, 2000, 5000 bp).</i>	E3161-01 E3161-02	100 loadings 500 loadings	€ 55.00 € 220.00
Perfect™ 100 bp DNA Ladder <i>Ready to Use DNA ladder for sizing DNA fragments from 100 to 2500 bp (100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 1500, 2000, 2500 bp).</i>	B E3134-01 E3134-02	50 µg 250 µg	€ 45.00 € 180.00
Perfect™ 100-1000 bp DNA Ladder <i>Ready to Use DNA ladder for sizing DNA fragments from 100 to 1000 bp (100, 200, 300, 400, 500, 600, 700, 800, 900, 1000 bp).</i>	E3141-01 E3141-02	50 µg 250 µg	€ 45.00 € 180.00
Perfect™ 1 kb DNA Ladder <i>Ready to Use DNA ladder for sizing DNA fragments from 500 to 10000 bp (500, 1000, 1500, 2000, 2500, 3000, 4000, 5000, 6000, 8000, 10000 bp).</i>	B E3130-01 E3130-02	100 µg 500 µg	€ 36.00 € 144.00
Perfect Plus™ 1 kb DNA Ladder <i>Ready to Use DNA ladder for sizing DNA fragments from 250 to 10000 bp (250, 500, 750, 1000, 1500, 2000, 2500, 3000, 4000, 5000, 6000, 8000, 10000 bp).</i>	B E3131-01 E3131-02	100 µg 500 µg	€ 39.00 € 156.00
FLUORESCENT-LABELED DNA MOLECULAR WEIGHT MARKERS			
GENEFLO™ 625 DNA ladder, ROX Labeled <i>Fluorescent-labeled DNA ladder for Use with ABI PRISM™ genetic analysers. Contains ROX dye.</i>	E3125-01 E3125-02	100 loadings 800 loadings	€ 50.00 € 250.00
GENEFLO™ 625 DNA ladder, TAMRA Labeled <i>Fluorescent-labeled DNA ladder for Use with ABI PRISM™ genetic analysers. Contains TAMRA dye.</i>	E3126-01 E3126-02	100 loadings 800 loadings	€ 50.00 € 250.00
5. MESOTHERMAL POLYMERASES			
DNA Polymerase Alpha <i>Human DNA Polymerase.</i>	E1075-01 E1075-02	50 U 200 U	€ 270.00 € 920.00
DNA Polymerase Gamma <i>Human DNA Polymerase.</i>	E1076-01 E1076-02	50 U 200 U	€ 270.00 € 920.00
DNA Polymerase Beta <i>Human DNA Polymerase.</i>	E1077-01 E1077-02	50 U 200 U	€ 210.00 € 720.00
Bst DNA Polymerase exo- <i>Large exonuclease free fragment of thermophilic Bst DNA Polymerase with strand displacement activity.</i>	E1078-01 E1078-02	1600 U 8000 U	€ 60.00 € 275.00
DNA Polymerase I <i>DNA-dependent DNA polymerase with inherent 3'→5' and 5'→3' exonuclease activity (E. coli).</i>	E1080-01 E1080-02	500 U 2 500 U	€ 45.00 € 180.00
Klenow Fragment <i>Large fragment of E. coli DNA Polymerase I (polymerase and proofreading 3'→5' exonuclease activities).</i>	E1091-01 E1091-02	200 U 1 000 U	€ 45.00 € 180.00
Klenow Fragment exo- <i>Large fragment of E. coli DNA Polymerase I (lacks 3'→5' exonuclease and 5'→3' exonuclease activities).</i>	E1092-01 E1092-02	200 U 1 000 U	€ 45.00 € 180.00


Vert!



PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
Poly (A) Polymerase <i>Poly (A) Polymerase is E. coli derived polymerase, polyadenylating 3' termini of RNA.</i>	E1240-01	40 U	€ 50.00
	E1240-02	200 U	€ 198.00
T4 DNA Polymerase <i>T4 DNA Polymerase exhibits very strong 3'->5' exonuclease activity.</i>	E1100-01	200 U	€ 54.00
	E1100-02	1 000 U	€ 216.00
T7 RNA Polymerase <i>Modified T7 RNA Polymerase with higher tolerance towards modified nucleotides.</i>	E1290-01	5 000 U	€ 68.00
	E1290-02	25 000 U	€ 225.00
Terminal Deoxynucleotidyl Transferase <i>Catalyzes the addition of deoxyribonucleotide monophosphate from triphosphate to the 3' hydroxyl terminus of DNA.</i>	E1390-01	500 U	€ 50.00
	E1390-02	2 500 U	€ 198.00
phi29 DNA Polymerase P <i>A highly processive recombinant polymerase with exceptional strand displacement activity.</i>	E1401-01	250 U	€ 45.00
	E1401-02	1 250 U	€ 180.00
MDA phi29 Kit P <i>Multiple displacement amplification (MDA) kit Using phi29 polymerase for whole-genome amplification experiments.</i>	E1402-01	25 reactions	€ 60.00
	E1402-02	125 reactions	€ 210.00

6. NUCLEASES

OMNI Nuclease <i>Nonspecific nuclease that completely degrades all forms of DNA and RNA.</i>	E1120-01	20 000 U	€ 30.00
	E1120-02	100 000 U	€ 135.00
DNase I (RNase-free) <i>Nonspecific DNase that degrades both double-stranded and single-stranded DNA endonucleolytically.</i>	E1345-01	1 000 U	€ 50.00
	E1345-02	5 000 U	€ 240.00
RNase I <i>Completely nonspecific RNase that hydrolyzes phosphodiester bond after all four bases.</i>	E1300-01	10 000 U	€ 63.00
	E1300-02	50 000 U	€ 252.00
RNase III <i>Cleaves long double-stranded RNA (dsRNA) into short dsRNAs (13-30 bases).</i>	E1340-01	200 U	€ 98.00
	E1340-02	1 000 U	€ 392.00
RNase A (DNase-free) <i>Endoribonuclease that attacks at the 3' phosphate of a pyrimidine nucleotide.</i>	E1350-01	10 mg (1 ml)	€ 23.00
	E1350-02	50 mg (5 ml)	€ 95.00
RNase H (E. coli) <i>Endonuclease that selectively cleaves the RNA of RNA-DNA hybrids.</i>	E1330-01	250 U	€ 49.50
	E1330-02	1 250 U	€ 198.00
Thermostable RNase H (Thermus thermophilus) <i>Endonuclease that selectively cleaves the RNA of RNA-DNA hybrids.</i>	E1325-01	500 U	€ 98.00
	E1325-02	2 500 U	€ 392.00
Exonuclease I B <i>Catalyzes the removal of nucleotides from single-stranded DNA in the 3' to 5' direction.</i>	E1150-01	4 000 U	€ 45.00
	E1150-02	20 000 U	€ 180.00
Exonuclease III <i>Exonuclease III is a 3'->5' exonuclease, releasing 5'-mononucleotides from the 3'-ends of DNA strands.</i>	E1140-01	25 000 U	€ 45.00
	E1140-02	125 000 U	€ 180.00
Exonuclease T5 <i>Catalyzes the removal of nucleotides from dsDNA and ssDNA in the 5' to 3' direction.</i>	E1130-01	1 000 U	€ 45.00
	E1130-02	5 000 U	€ 180.00
Lambda Exonuclease <i>Double-stranded specific DNase, digests one DNA strand starting from 5' terminus.</i>	E1180-01	1 000 U	€ 45.00
	E1180-02	5 000 U	€ 180.00
Mung Bean Nuclease <i>Single-stranded specific DNA and RNA endonuclease.</i>	E1190-01	2 000 U	€ 45.00
	E1190-02	10 000 U	€ 180.00
S1 Nuclease <i>Non-specific endonuclease, active primarily on single-stranded DNA and RNA.</i>	E1335-01	10 000 U	€ 36.00
	E1335-02	50 000 U	€ 141.00
T7 Endonuclease I <i>Recognizes and cleaves non-perfectly matched DNA, cruciform DNA structures, Holliday structures or junctions, heteroduplex DNA and more slowly, nicked ds-DNA.</i>	E1125-01	250 U	€ 65.00
	E1125-02	1 250 U	€ 290.00



new product



bestseller



unique



promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
7. LIGASES & CLONING KITS			
DNA Ligase (<i>E. coli</i>)	E1065-01	500 U	€ 113.00
	E1065-02	2 500 U	€ 450.00
T4 DNA Ligase (<i>E. coli</i> infected with bacteriophage T4)	E1060-01	20 000 U	€ 54.00
	E1060-02	100 000 U	€ 216.00
Taq DNA Ligase (<i>Thermus aquaticus</i>)	E1070-01	1 000 U	€ 36.00
	E1070-02	5 000 U	€ 144.00
LigON Kit	E1050-01	10 reactions	€ 100.00
	E1050-02	50 reactions	€ 400.00
SpeedyLink DNA Ligation Kit	E1055-01	30 reactions	€ 75.00
	E1055-02	150 reactions	€ 300.00
<i>Kit enables ligation of cohesive end or blunt end DNA fragments in 5 minutes at room temperature.</i>			

8. SPECIAL PURPOSE ENZYMES & REAGENTS


DNA Modification			
Alkaline Phosphatase (<i>Calf Intestinal</i>)	E1025-01	1 000 U	€ 25.00
	E1025-02	5 000 U	€ 105.00
<i>Calf intestine derived phosphatase which catalyzes the the removal of 5'-phosphate from DNA.</i>			
Polar - BAP Thermosensitive Alkaline Phosphatase	E1027-01	1 000 U	€ 48.00
	E1027-02	5 000 U	€ 235.00
<i>Catalyzes the removal of 5'-phosphate from DNA and RNA. Heat inactivation (5 minutes at 70°C).</i>			
Fast Polar - BAP Thermosensitive Alkaline Phosphatase	E1028-01	1 000 U	€ 58.00
	E1028-02	5 000 U	€ 280.00
<i>Catalyzes the removal of 5'-phosphate from DNA and RNA. Heat inactivation (5 minutes at 70°C), short reaction time.</i>			
Exo-BAP Kit*	E1011-01	100 reactions	€ 35.00
	E1011-02	500 reactions	€ 170.00
<i>* Please check the availability for your country.</i>			
<i>Allows quick and effective removing of any Unconsumed dNTPs and primers prior sequencing reaction (contains enzymes in separate vials).</i>			
Exo-BAP Mix*	E1012-01	100 reactions	€ 35.00
	E1012-02	500 reactions	€ 170.00
<i>* Please check the availability for your country.</i>			
<i>Allows quick and effective removing of any Unconsumed dNTPs and primers prior sequencing reaction (contains optimized mixture of enzymes).</i>			
T4 Polynucleotide Kinase	E1261-01	1 000 U	€ 68.00
	E1261-02	5 000 U	€ 272.00
<i>Catalyzes the phosphorylation of 5' hydroxyl termini of ds and ss DNA or RNA.</i>			
Thermophilic Pyrophosphatase (<i>Thermus aquaticus</i>)	E1267-01	50 U	€ 68.00
	E1267-02	250 U	€ 272.00
<i>Thermostable pyrophosphatase, hydrolyzing inorganic pyrophosphate into phosphate.</i>			
Rec A Protein (<i>E. coli</i>)	E1370-01	0.2 mg	€ 62.00
	E1370-02	1 mg	€ 245.00
<i>Recombined protein that enhances interaction and exchange of homologous DNA strands.</i>			
Endocyclase N15 (<i>Protelomerase N15</i>)	E4410-01	50 U	€ 18.00
	E4410-02	250 U	€ 75.00
<i>Catalyzes formation of covalent links between complementary DNA 5'- and 3'- ends (hairpin).</i>			

DNA Replication

SV40 Large Antigen	E5800-01	20 µg	€ 170.00
	E5800-02	100 µg	€ 650.00
<i>The major regulatory protein of SV40 virus.</i>			

COT DNA

Human COT DNA	E4110-01	0.5 mg	€ 95.00
	E4110-02	2 mg	€ 360.00
<i>Human placenta COT DNA optimized for a size range of 50 to 300 bp.</i>			
Uncut Human Placental DNA	E4111-01	0.5 mg	€ 95.00
	E4111-02	2 mg	€ 360.00



new product



bestseller




unique



promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
DNA - precipitation & visualisation			
Vivid Violet DNA/ RNA Co-Precipitate B	E4502-01	250 µl	€ 45.00
SimplySafe™ B <i>New and safe nucleic acid stain, an alternative to the traditional ethidium bromide stain.</i>	E4600-01 E4600-02	1 ml 10 ml	€ 48.00 € 405.00
Enzyme stabilizing agents			
Bovine Serum Albumin <i>BSA acetylated to inactivate contaminating nucleases. Used as a supplement in biochemical reactions.</i>	E4020-01 E4020-02	75 mg 375 mg	€ 45.00 € 180.00
Agarases			
Agarase <i>Allows for a simple, quantitative recovery of intact nucleic acids from agarose gels.</i>	E4800-01 E4800-02	100 U 500 U	€ 50.00 € 180.00
9. WATER, BUFFERS & SOLUTIONS (Molecular Biology Grade)			
RNase-free water <i>DEPC-treated and nuclease-free water for RNA work.</i>	E0210-01 E0210-02 E0210-03	5 x 1.5 ml 5 x 5.0 ml 500 ml	€ 18.00 € 31.00 € 39.00
PCR water <i>Ultrapure, sterile, deionized, double distilled and prealiquoted water.</i>	E0211-01 E0211-02	5 x 1.5 ml 10 x 1.5 ml	€ 25.00 € 45.00
Nuclease-free water, not DEPC treated <i>For DNA targeted reactions such as PCR, ligation and restriction digestions.</i>	E0212-01 E0212-02	100 ml 1 l	€ 18.00 € 48.00
50 x TAE	E0220-01 E0220-02	500 ml 1 l	€ 30.00 € 50.00
5 x TBE	E0230-01 E0230-02	1 l 5 x 1 l	€ 20.00 € 90.00
0.5 M EDTA pH 8.0	E0240-01 E0240-02 E0240-03	100 ml 500 ml 1 l	€ 16.00 € 25.00 € 45.00
1 x TE pH 8.0	E0250-01 E0250-02 E0250-03	500 ml 1 l 5 x 1 l	€ 16.00 € 22.00 € 100.00
1M Tris-HCl pH 6.8	E0270-01	1 l	€ 35.00
1M Tris-HCl pH 7.0	E0271-01	1 l	€ 35.00
1M Tris-HCl pH 7.5	E0272-01	1 l	€ 35.00
1M Tris-HCl pH 8.0	E0273-01	1 l	€ 35.00
1M Tris-HCl pH 8.8	E0274-01	1 l	€ 35.00
5 x TG	E0254-01 E0254-02	500 ml 1 l	€ 30.00 € 55.00
5 x TGS	E0255-01 E0255-02	500 ml 1 l	€ 30.00 € 55.00
1 x PBS	E0281-01 E0281-02	500 ml 1 l	€ 12.00 € 20.00
1 x PBS (without Mg and Ca ions)	E0282-01 E0282-02	500 ml 1 l	€ 10.00 € 18.00
6 x Loading Buffer BLUE <i>For loading a DNA sample on agarose and polyacrylamide gels.</i>	E0260-01	5 x 1 ml	€ 25.00

Verte!




new product



bestseller



unique



promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
6 x Loading Buffer RED <i>For loading a DNA sample on agarose and polyacrylamide gels.</i>	E0261-01	5 x 1 ml	€ 25.00
6 x Loading Buffer GREEN <i>For loading a DNA sample on agarose and polyacrylamide gels.</i>	E0262-01	5 x 1 ml	€ 25.00
6 x Loading Buffer YELLOW <i>For loading a DNA sample on agarose and polyacrylamide gels.</i>	E0263-01	5 x 1 ml	€ 25.00
6 x Loading Buffer TriDye <i>For loading a DNA sample on agarose and polyacrylamide gels.</i>	E0265-01	5 x 1 ml	€ 25.00
Extraction Buffer (CTAB based) <i>For quick and efficient lysis of plant cells.</i>	E0290-01	1 l	€ 50.00
RIPA Lysis Buffer <i>RIPA (Radioimmunoprecipitation Assay) Buffer for general lysis of cells.</i>	E0295-01 E0295-02	50 ml 200 ml	€ 38.00 € 170.00
BlueStain Sensitive* <i>* Please check the availability for your country Fast and safe staining of proteins after SDS-PAGE. Minimal sensitivity - 35 ng of protein.</i>	E0298-01	1 l	€ 44.00
BlueStain Sensitive Plus* <i>* Please check the availability for your country. Fast and safe staining of proteins after SDS-PAGE. Minimal sensitivity - 20 ng of protein.</i>	E0297-01	1 l	€ 42.00
4 x Protein Loading Buffer <i>SDS-PAGE sample buffer, without β-mercaptoethanol.</i>	E0268-01	10 x 1 ml	€ 25.00
4 x Protein Loading Buffer Plus <i>SDS-PAGE sample buffer, with β-mercaptoethanol.</i>	E0269-01	10 x 1 ml	€ 25.00

10. XELEX KITS

Aptamers allow distinct molecules to bind to target structures with high specificity and affinity by non-covalent molecular interactions (with respect to certain features aptamers resemble antibodies). Aptamer selection against a broad range of target structures is possible, including proteins, polysaccharides, various large and small organic molecules (drug compounds) and even small metal ions.




XELEX Kits for aptamer in vitro evolution and selection.

XELEX DNA CORE Kit (Combo)	E3650-01	50 sel. / 25 anal.	€ 245.00
XELEX DNA CORE Kit (Selection Unit)	E3651-01	50 selections	€ 155.00
XELEX DNA CORE Kit (Analysis Unit)	E3652-01	25 analysis	€ 100.00
XELEX RNA CORE Kit	E3660-01	50 selections	€ 365.00

11. PROTEASES, ENZYMES & REAGENTS FOR PROTEOMICS

Proteinase K <i>Solution 20 mg/ml.</i>	E4350-01 E4350-02	1 ml 5 x 1 ml	€ 20.00 € 90.00
Sortase A <i>Protein-protein ligase/ transpeptidase.</i>	E4400-01	20 μ g	€ 75.00
TEV Protease <i>A cysteine protease that specifically recognizes and cleaves a linear epitope.</i>	E4310-01 E4310-02	1 000 U 10 000 U	€ 60.00 € 540.00
Lysozyme <i>Improves bacterial cell lysis (10 mg/ml).</i>	E4320-01 E4320-02	1 ml 5 x 1 ml	€ 25.00 € 120.00



PRODUCT NAME & DESCRIPTION

CAT.NUMBER

PACKAGE

PRICE

12. TISSUE CULTURE MEDIA COMPONENTS & REAGENTS

Fetal Bovine Serum*Origin: Brazil, certified for EU (not heat inactivated).*

E5050-02

500 ml

call for price

Fetal Bovine Serum HI*Origin: Brazil, certified for EU (heat inactivated).*

E5051-02

500 ml

call for price

13. AGAROSES (Molecular Biology Grade)

Agarose*Ideal for everyday analysis of nucleic acids by gel electrophoresis or blotting.*

E0301-100

100 g

€ 40.00

E0301-500

500 g

€ 150.00

High ReResolution Agarose*Separates small DNA and RNA fragments between 20 and 800 bp in length.*

E0302-50

50 g

€ 120.00

Low Melting Agarose*Recommended for convenient recovery of nucleic acids after electrophoresis.*

E0303-50

50 g

€ 115.00

Agarose tablets*Agarose is supplied in convenient tablet form. Each agarose tablet consists 0.5 g agarose.*

E0305-01

50 tablets

€ 30.00

E0305-02

200 tablets

€ 110.00

14. RESTRICTION ENDONUCLEASES

Enzyme	Prototype	Reaction temperature	Recognition sequence	Catalog number	Package	Price
AccI	AccI	37°C	GT / MKAC	E2004-01	500 U	€ 50.00
				E2004-02	2 500 U	€ 200.00
AcvI	PmaCI	37°C	CAC / GTG	E2005-01	1 500 U	€ 36.00
				E2005-02	7 500 U	€ 143.00
AluI	AluI	37°C	AG / CT	E2010-01	500 U	€ 33.00
				E2010-02	2 500 U	€ 130.00
ApaI	ApaI	25°C	GGCC / C	E2015-01	5 000 U	€ 45.00
				E2015-02	25 000 U	€ 180.00
AvaI	AvaI	37°C	C / YCGRG	E2006-01	2 000 U	€ 36.00
				E2006-02	10 000 U	€ 143.00
Avall	Avall	37°C	G / GWCC	E2007-01	1 000 U	€ 39.00
				E2007-02	5 000 U	€ 156.00
Ball	Ball	37°C	TGG / CCA	E2008-01	250 U	€ 57.00
				E2008-02	1 250 U	€ 228.00
BamHI B	BamHI	37°C	G / GATCC	E2050-01	10 000 U	€ 36.00
				E2050-02	50 000 U	€ 143.00
BanII	HgiIII	37°C	GRGCY / C	E2060-01	1 000 U	€ 41.00
				E2060-02	5 000 U	€ 162.00
BglI	BglI	37°C	GCCNNNN / NGGC	E2080-01	2 000 U	€ 36.00
				E2080-02	10 000 U	€ 143.00
BglII	BglII	37°C	A / GATCT	E2090-01	2 000 U	€ 36.00
				E2090-02	10 000 U	€ 143.00
BsiHKCI	AvaI	65°C	C / YCGRG	E2109-01	10 000 U	€ 41.00
				E2109-02	50 000 U	€ 162.00
BssHII	BsePI	50°C	G / CGCGC	E2115-01	200 U	€ 41.00
				E2115-02	1 000 U	€ 162.00
BstXI	BstXI	50°C	CCANNNNN / NTGG	E2118-01	1 000 U	€ 41.00
				E2118-02	5 000 U	€ 162.00
Clal	Clal	37°C	AT / CGAT	E2110-01	1 000 U	€ 46.00
				E2110-02	5 000 U	€ 184.00
CviJI	CviJI	37°C	RG / CY	E2125-01	100 U	€ 113.00
				E2125-02	400 U	€ 360.00
CviJI* U	CviJI*	37°C	G / C except YG / CR	E2126-01	100 U	€ 149.00
				E2126-02	400 U	€ 475.00
Ddel	Ddel	37°C	C / TNAG	E2140-01	500 U	€ 50.00
				E2140-02	2 500 U	€ 200.00

Vertel!

N new productB bestsellerU uniqueP promotion

PRODUCT NAME & DESCRIPTION				CAT.NUMBER	PACKAGE	PRICE
Dpnl	Dpnl	37°C	GA / TC	E2135-01	500 U	€ 35.00
				E2135-02	2 500 U	€ 135.00
Dral	AhaIII	37°C	TTT / AAA	E2145-01	2 000 U	€ 33.00
				E2145-02	10 000 U	€ 130.00
EcoRI B	EcoRI	37°C	G / AATTC	E2150-01	10 000 U	€ 33.00
				E2150-02	50 000 U	€ 130.00
EcoRV	EcoRV	37°C	GAT / ATC	E2160-01	2 000 U	€ 33.00
				E2160-02	10 000 U	€ 130.00
FokI	FokI	37°C	GGATG (9/13)	E2170-01	500 U	€ 36.00
				E2170-02	2 500 U	€ 143.00
HaeIII	HaeIII	37°C	GG / CC	E2180-01	3 000 U	€ 33.00
				E2180-02	15 000 U	€ 130.00
HhaI	HhaI	37°C	GCG / C	E2190-01	2 000 U	€ 50.00
				E2190-02	10 000 U	€ 200.00
HincII	HindII	37°C	GTY / RAC	E2200-01	1 000 U	€ 36.00
				E2200-02	5 000 U	€ 143.00
HindIII B	HindIII	37°C	A / AGCTT	E2220-01	10 000 U	€ 33.00
				E2220-02	50 000 U	€ 130.00
Hinfl	Hinfl	37°C	G / ANTC	E2240-01	5 000 U	€ 36.00
				E2240-02	25 000 U	€ 143.00
HpaI	HpaI	37°C	GTT / AAC	E2250-01	500 U	€ 33.00
				E2250-02	2 500 U	€ 130.00
HpaII	HpaII	37°C	C / CGG	E2260-01	2 000 U	€ 33.00
				E2260-02	10 000 U	€ 130.00
KpnI	KpnI	37°C	GGTAC / C	E2280-01	3 000 U	€ 36.00
				E2280-02	15 000 U	€ 143.00
MboI	MboI	37°C	/ GATC	E2283-01	200 U	€ 36.00
				E2283-02	1 000 U	€ 143.00
MboII	MboII	37°C	GAAGA (8/7)	E2284-01	250 U	€ 36.00
				E2284-02	1 250 U	€ 143.00
MluI N	MluI	37°C	A / CGCGT	E2288-01	1 000 U	€ 36.00
				E2288-02	5 000 U	€ 143.00
MnII	MnII	37°C	CCTC (7/6)	E2289-01	250 U	€ 36.00
				E2289-02	1 250 U	€ 143.00
MspI	HpaII	37°C	C / CGG	E2290-01	5 000 U	€ 36.00
				E2290-02	25 000 U	€ 143.00
NarI	NarI	37°C	GG / CGCC	E2291-01	200 U	€ 36.00
				E2291-02	1 000 U	€ 143.00
NcoI B	NcoI	37°C	C / CATGG	E2292-01	1 000 U	€ 36.00
				E2292-02	5 000 U	€ 143.00
NdeI	NdeI	37°C	CA / TATG	E2293-01	500 U	€ 25.00
				E2293-02	2 500 U	€ 100.00
NheI	NheI	37°C	G / CTAGT	E2294-01	500 U	€ 25.00
				E2294-02	2 500 U	€ 100.00
NotI	NotI	37°C	GC / GGCCGC	E2296-01	500 U	€ 50.00
				E2296-02	2 500 U	€ 200.00
NruI	NruI	37°C	TCG / CGA	E2297-01	1 000 U	€ 36.00
				E2297-02	5 000 U	€ 143.00
NsiI	NsiI	37°C	ATGCA / T	E2298-01	1 000 U	€ 50.00
				E2298-02	5 000 U	€ 200.00
PinAI	AgeI	37°C	A / CCGGT	E2299-01	200 U	€ 55.00
				E2299-02	1 000 U	€ 217.00
PstI	PstI	37°C	CTGCA / G	E2300-01	10 000 U	€ 41.00
				E2300-02	50 000 U	€ 162.00
PvuI	PvuI	37°C	CGAT / CG	E2320-01	200 U	€ 33.00
				E2320-02	1 000 U	€ 130.00
PvuII	PvuII	37°C	CAG / CTG	E2331-01	5 000 U	€ 41.00
				E2331-02	25 000 U	€ 162.00
RsaI	RsaI	37°C	GT / AC	E2340-01	1 000 U	€ 36.00
				E2340-02	5 000 U	€ 143.00
RsrII	RsrII	37°C	CG / GWCCG	E2342-01	100 U	€ 36.00
				E2342-02	500 U	€ 143.00
SacI	SacI	37°C	GAGCT / C	E2350-01	2 000 U	€ 36.00
				E2350-02	10 000 U	€ 143.00

Verte!

N new product

B bestseller

U unique

P promotion

PRODUCT NAME & DESCRIPTION				CAT.NUMBER	PACKAGE	PRICE	
SaclI	SaclI	37°C	CCGC / GG	E2352-01	2 000 U	€ 36.00	
				E2352-02	10 000 U	€ 143.00	
Sall	Sall	37°C	G / TCGAC	reduced star activity	E2371-01	2 000 U	€ 41.00
					E2371-02	10 000 U	€ 162.00
Sau3AI	Mbol	37°C	/ GATC	E2375-01	200 U	€ 41.00	
				E2375-02	1 000 U	€ 162.00	
Scal	Scal	37°C	AGT / ACT	E2378-01	1 000 U	€ 36.00	
				E2378-02	5 000 U	€ 143.00	
Scal	Scal	37°C	AGT / ACT	reduced star activity	E2379-01	1 000 U	€ 41.00
					E2379-02	5 000 U	€ 162.00
Sfil	Sfil	50°C	GGCCNNNN / NGGCC	E2380-01	2 000 U	€ 41.00	
				E2380-02	10 000 U	€ 162.00	
SmaI	SmaI	25°C	CCC / GGG	E2390-01	2 000 U	€ 36.00	
				E2390-02	10 000 U	€ 143.00	
SpeI	SpeI	37°C	A / CTAGT	E2398-01	500 U	€ 41.00	
				E2398-02	2 500 U	€ 162.00	
SphI	SphI	37°C	GCATG / C	E2400-01	250 U	€ 36.00	
				E2400-02	1 250 U	€ 143.00	
SspI	SspI	37°C	AAT / ATT	E2402-01	1 000 U	€ 41.00	
				E2402-02	5 000 U	€ 162.00	
StuI	StuI	37°C	AGG / CCT	E2405-01	1 000 U	€ 36.00	
				E2405-02	5 000 U	€ 143.00	
TaqI	TaqI	65°C	T / CGA	E2410-01	4 000 U	€ 36.00	
				E2410-02	20 000 U	€ 143.00	
TaqII	U TaqII	65°C	GACCGA (11/9)	E2411-01	100 U	€ 81.00	
				E2411-02	500 U	€ 324.00	
TspDTI	U TspDTI	70°C	ATGAA (11/9)	E2502-01	50 U	€ 81.00	
				E2502-02	250 U	€ 324.00	
TspGWI	U TspGWI	70°C	ACGGA (11/9)	E2501-01	50 U	€ 81.00	
				E2501-02	250 U	€ 324.00	
Tth111I	Tth111I	65°C	GACN / NNGTC	E2420-01	1 000 U	€ 45.00	
				E2420-02	5 000 U	€ 180.00	
XbaI	XbaI	37°C	T / CTAGA	E2430-01	3 000 U	€ 36.00	
				E2430-02	15 000 U	€ 143.00	
XhoI	XhoI	37°C	C / TCGAG	E2440-01	4 000 U	€ 36.00	
				E2440-02	20 000 U	€ 143.00	

15. RNA PURIFICATION SYSTEMS

GeneMATRIX series

GeneMATRIX ENVIRONMENTAL DNA & RNA Purification Kit <i>Kit for isolation and purification genomic DNA and total RNA. For stool, soil and environmental samples.</i>	E3572-01	25 preps	€ 58.00
	E3572-02	100 preps	€ 215.00
GeneMATRIX VIRAL RNA/DNA Purification Kit <i>Kit for simultaneous purification of viral RNA / DNA from plasma, serum and cell-free body fluids.</i>	E3592-01	25 preps	€ 58.00
	E3592-02	100 preps	€ 215.00
GeneMATRIX FFPE RNA Purification Kit <i>Kit for isolation of RNA from formalin-fixed, paraffin-embedded samples.</i>	E3593-01	25 preps	€ 68.00
	E3593-02	100 preps	€ 250.00
GeneMATRIX UNIVERSAL BLOOD RNA Purification Kit <i>Universal Kit for isolation of total RNA from human or animal blood sample, fresh or frozen with stabilizers.</i>	E3594-01	25 preps	€ 132.00
	E3596-01	25 preps	€ 58.00
GeneMATRIX HUMAN BLOOD RNA Purification Kit <i>Kit for isolation of total RNA from fresh human blood from leukocytes.</i>	E3597-01	25 preps	€ 68.00
	E3597-02	100 preps	€ 250.00
GeneMATRIX UNIVERSAL DNA/RNA/Protein Purification Kit <i>Kit for isolation of total DNA/RNA/Protein from the same biological sample (bacteria, tissue, plant, yeast, cell culture).</i>	E3598-01	25 preps	€ 52.00
	E3598-02	100 preps	€ 195.00

Vertel!

N new product

B bestseller

U unique

P promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
GeneMATRIX UNIVERSAL RNA/miRNA Purification Kit <i>Universal Kit for isolation of miRNA, total RNA or miRNA and total RNA together, from the tissue, plant and cell culture, as well as for purification of miRNA after enzymatic reactions.</i>	E3599-01 E3599-02	25 preps 100 preps	€ 58.00 € 215.00
GeneMATRIX Tick RNA/DNA Purification Kit <i>Kit for simultaneous isolation of viral and bacterial RNA and DNA that live in small arachnids or insects.</i>	E3590-01	25 preps	€ 68.00
Universal RNA Clean-Up Kit <i>Kit removes DNA and other contaminants from RNA samples isolated from different sources (animal and plant tissues, bacteria, yeast, blood etc.).</i>	E3589-01	25 preps	€ 52.00

GeneMAGNET series


GeneMAGNET RNA Purification Kit <i>Kit for isolation of total RNA from various sample (animal and plant tissue, bacteria, yeast, cell culture), with magnetic beads.</i>	E3401-01	96 preps	€ 240.00
GeneMAGNET Viral RNA/DNA Purification Kit <i>Kit for simultaneous purification of viral RNA / DNA from plasma, serum and cell-free body fluids, with magnetic beads.</i>	E3402-01	96 preps	€ 250.00

16. DNA PURIFICATION SYSTEMS

GeneMATRIX series

GeneMATRIX PLASMID MINIPREP DNA Purification Ki <i>Kit for isolation of high-purity plasmid DNA (1.5-4 ml bacterial culture).</i>	E3500-01 E3500-02	50 preps 150 preps	€ 49.00 € 138.00
GeneMATRIX BIO-TRACE DNA Purification Kit B <i>Kit for isolation of DNA from biological traces from various samples for clinical and forensic analysis.</i>	E3510-01 E3510-02	25 preps 100 preps	€ 54.00 € 200.00
GeneMATRIX SHORT DNA CLEAN - UP Purification Kit <i>Kit for purification of short single-stranded and double-stranded DNA fragments after enzymatic reactions.</i>	E3515-01 E3515-02	25 preps 100 preps	€ 39.00 € 140.00
GeneMATRIX PCR/ DNA CLEAN - UP Purification Kit <i>Kit for purification of PCR products / DNA after enzymatic reactions.</i>	E3520-01 E3520-02	50 preps 150 preps	€ 42.00 € 115.00
GeneMATRIX FOOD-EXTRACT DNA Purification Kit <i>Kit for isolation of DNA from fresh and processed food of plant or animal origin.</i>	E3525-01 E3525-02	25 preps 100 preps	€ 54.00 € 200.00
GeneMATRIX SWAB-EXTRACT DNA Purification Kit <i>Kit for isolation of DNA from swabs for clinical and forensic analysis.</i>	E3530-01 E3530-02	25 preps 100 preps	€ 39.00 € 145.00
GeneMATRIX AGROBACTERIUM PLASMID MINIPREP DNA Purification Kit <i>Kit for isolation of DNA from Agrobacterium tumefaciens (0.5-2 ml bacterial culture).</i>	E3535-01 E3535-02	50 preps 150 preps	€ 64.00 € 180.00
GeneMATRIX AGAROSE - OUT DNA Purification Kit B <i>Universal Kit for purification of DNA from agarose gels.</i>	E3540-01 E3540-02	50 preps 150 preps	€ 50.00 € 140.00
GeneMATRIX BASIC DNA Purification Kit B <i>Universal Kit for purification of PCR products / DNA after enzymatic reactions, isolation of DNA from agarose gels and purification of plasmid DNA from bacteria.</i>	E3545-01 E3545-02	50 preps 150 preps	€ 52.00 € 154.00
GeneMATRIX TISSUE DNA Purification Kit <i>Universal Kit for isolation of DNA from human and animal tissues, cultured cells, rodents tails, hair and insects.</i>	E3550-01 E3550-02	50 preps 150 preps	€ 80.00 € 225.00
GeneMATRIX TISSUE & BACTERIAL DNA Purification Kit <i>Kit for isolation of total DNA from human and animal tissues and bacteria.</i>	E3551-01 E3551-02	50 preps 150 preps	€ 82.00 € 230.00
GeneMATRIX CELL CULTURE DNA Purification Kit <i>Kit for isolation of DNA from human and animal cell cultures.</i>	E3555-01 E3555-02	50 preps 150 preps	€ 80.00 € 225.00

Verte!



PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
GeneMATRIX BONE DNA Purification Kit <i>Kit for isolation of DNA from animal or human bones.</i>	E3560-01 E3560-02	25 preps 50 preps	€ 73.00 € 135.00
GeneMATRIX QUICK BLOOD DNA Purification Kit <i>Kit for quick isolation of DNA from fresh or frozen blood.</i>	B E3565-01 E3565-02	50 preps 150 preps	€ 86.00 € 235.00
GeneMATRIX SOIL DNA Purification Kit <i>Kit for isolation of DNA from soil.</i>	E3570-01 E3570-02	50 preps 100 preps	€ 104.00 € 200.00
GeneMATRIX STOOL DNA Purification Kit <i>Kit for isolation of DNA from stool.</i>	E3575-01 E3575-02	50 preps 100 preps	€ 104.00 € 200.00
GeneMATRIX BACTERIAL & YEAST GENOMIC DNA Purification Kit <i>Universal Kit for isolation of total DNA from bacteria Gram+, Gram- and yeast.</i>	E3580-01 E3580-02	50 preps 150 preps	€ 80.00 € 225.00
GeneMATRIX GRAM PLUS & YEAST GENOMIC DNA Purification Kit <i>Universal Kit for isolation of total DNA from bacteria Gram+ and yeast.</i>	E3585-01 E3585-02	25 preps 100 preps	€ 49.00 € 190.00
GeneMATRIX PLANT & FUNGI DNA Purification Kit <i>Universal Kit for isolation of total DNA from plants, fungi and lichens.</i>	E3595-01 E3595-02	50 preps 150 preps	€ 92.00 € 255.00
DNA/RNA Extracol Kit <i>Complete Kit containing spin columns and phenol based chemicals for the isolation of total DNA and RNA from diverse samples.</i>	E3750-01 E3750-02	25 preps 100 preps	€ 50.00 € 175.00
Universal Kit for DNA isolation with GeDI Solution <i>Complete Kit containing spin columns and ready-to-use GeDI reagent for the isolation of total DNA</i>	E3765-01 E3765-02	25 preps 100 preps	€ 50.00 € 175.00
AnEx Plasmid Midiprep DNA purification Kit <i>Kit for isolation of plasmid DNA in midi format based on ion exchange chromatography.</i>	E3790-01 E3790-02	4 preps 8 preps	€ 25.00 € 48.00
GeneMAGNET series			
GeneMAGNET PCR/ DNA CLEAN - UP Purification Kit <i>Kit for purification of PCR products / DNA after enzymatic reactions, with magnetic beads.</i>	E3420-01	96 preps	€ 115.00
GeneMAGNET Human and Animal Tissue DNA Purification Kit <i>Kit for isolation of DNA from human and animal tissues, with magnetic beads.</i>	N E3425-01	96 preps	€ 240.00
GeneMAGNET Food DNA Purification Kit <i>Kit for isolation of DNA from fresh and processed food of plant or animal origin, with magnetic beads.</i>	N E3427-01	96 preps	€ 240.00
17. EMULSION KITS			
Micellula DNA Emulsion & Purification Kit <i>Kit for creation of emulsions and subsequent DNA purification. For in vitro compartmentalization (IVC) and for emulsion PCR (ePCR).</i>	U E3600-01 E3600-02	50 preps 150 preps	€ 75.00 € 210.00

PRODUCT NAME & DESCRIPTION

CAT.NUMBER

PACKAGE

PRICE

18. CHEMICALS, BUFFERS & ACCESSORIES FOR DNA/RNA PURIFICATION SYSTEMS

Buffer from DNA/RNA Purification Kit <i>Not applicable for AnEx Plasmid Midiprep DNA Purification Kit, DNA/RNA Extracol Kit and LyseBlood buffer.</i> <i>While ordering please indicate the catalog number of the DNA/RNA Purification System.</i>	E0309-01	25/50 preps	€ 12.00
	E0309-02	100/150 preps	€ 25.00
RNA Storage N <i>For precipitate dissolving and long-term storage of RNA.</i>	E0308-01	10x1,5 ml	€ 55.00
Lyse CT buffer <i>Solely for E3595 system. Additional lysing Solution for plant tissues rich in polysaccharides.</i>	E0324-01	50 preps	€ 11.00
	E0324-02	150 preps	€ 25.00
Lyse Buffer PVP N <i>Solely for E3597 and E3598 systems. Additional lysing Solution for plant tissues rich in polysaccharides and phenolic compounds.</i>	E0291-01	10 ml	€ 35.00
Lyse RBC buffer 1x <i>Redy for Use lysing Solution for red blood cells.</i>	E0326-01	100 ml	€ 10.00
	E0326-02	250 ml	€ 20.00
fix RNA* B <i>* Please check the availability for your country.</i> <i>RNA stabilization reagent.</i>	E0280-01	100 ml	€ 48.00
	E0280-02	200 ml	€ 90.00
	E0280-03	500 ml	€ 215.00
RNA Extracol <i>Ready-to-use reagent, phenol based, for the isolation of total RNA from samples of diverse origins.</i>	E3700-01	25 ml	€ 15.00
	E3700-02	100 ml	€ 55.00
AFR01 <i>Antifoaming reagent for EURx lysis buffers.</i>	E0328-01	0,2 ml	€ 10.00
	E0328-02	1 ml	€ 45.00
Verte!			
BeadTubeDry <i>Tubes containing glass spheres in two sizes for homogenization samples of diverse origin, especially bacteria and yeasts.</i>	E0358-01	25 preps	€ 25.00
	E0358-02	100 preps	€ 80.00
Tissue Grinding Tool <i>Combination of microtube, beads and pestle for homogenization of small samples of tissue, bacteria, yeast or cell cultures. One can choose from two types of beads.</i> <i>While ordering please indicate the a/b letter next to a catalog number.</i>	E0359-01	10 preps	€ 15.00
	E0359-02	50 preps	€ 65.00
Universal Solution for DNA isolation (GeDI) <i>Ready-to-use reagent, mixture of detergents and salts for DNA isolation from samples of diverse origins.</i>	E3760-01	25 ml	€ 25.00
	E3760-02	100 ml	€ 80.00
Tissue mini grinder <i>Ready-to-use tissue mini grinder.</i>	E0357-01	1 piece	€ 50.00

N new product

B bestseller

U unique

P promotion

PRODUCT NAME & DESCRIPTION	CAT.NUMBER	PACKAGE	PRICE
19. COLUMNS FOR DNA/RNA PURIFICATION SYSTEMS			
Columns for DNA Purification Columns for selected DNA Purification Systems. <i>While ordering please indicate the catalog number of the DNA Purification System.</i>	E0330-01 E0330-02	25 columns 100 columns	€ 20.00 € 70.00
Columns for DNA Purification N Columns for selected DNA Purification Systems, individually packed. <i>While ordering please indicate the catalog number of the DNA Purification System.</i>	E0331-01	100 columns	€ 80.00
Columns for RNA Purification Columns set (binding and homogenization columns) for selected RNA Purification Systems. <i>While ordering please indicate the catalog number of the RNA Purification System.</i>	E0340-01 E0340-02	2x25 columns 2x100 columns	€ 30.00 € 110.00
Columns for RNA Purification N Columns set (binding and homogenization columns) for selected RNA Purification Systems, individually packed. <i>While ordering please indicate the catalog number of the RNA Purification System.</i>	E0341-01	2x50 columns	€ 80.00
Collecting tubes without caps Collecting tubes (2.0 ml) without caps, ready to Use for all purification steps except elution. <i>While ordering please indicate the catalog number of the DNA/RNA Purification System.</i>	E0350-01 E0350-02	25 tubes 100 tubes	€ 20.00 € 50.00
Elution tubes Elution tubes (1.5 ml) ready for elution of DNA and RNA. Capped tubes to ensure no loss of sample.	E0355-01 E0355-02	25 tubes 100 tubes	€ 20.00 € 50.00



EURx Sp. z o.o.
tel. +48 58 5240697
fax +48 58 3417423
www.eurx.com.pl
sales@eurx.com.pl