



PCR Water

(PCR Grade)

Cat. No.	size
E0211-01	5 x 1.5 ml
E0211-02	10 x 1.5 ml

Storage Conditions:

Store at -20°C (long term) or at 2-8°C (short term).

Ultrapure, sterile, deionized, destiled and prealiquoted water.

Suitable for use in:

- all PCRs and RT-PCRs,
- all conventional and qPCRs with any instrument,
- cDNA synthesis,
- sequencing / cycle sequencing.

Quality Control:

All lots are monitored for contaminating 3'-exonuclease, 5'-exonuclease as well as for nonspecific single- and double-stranded DNase activities.

All lots are PCR-assayed for absence of contaminating DNA from bacterial, fungi and human sources.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Sp. z o.o. ul. Przyrodników 3, 80-297 Gdańsk

NIP 957-07-05-191, KRS 0000202039, www.eurx.com.pl, orders@eurx.com.pl

tel. +48 58 524 06 97, fax +48 58 341 74 23