



50 x TAE

Molecular Biology Grade

Cat No.	size
E0220-01	500 ml
E0220-02	1 l

Storage:

Store at room temperature.

Description:

50 x TAE Buffer is a sterile-filtered solution of 2 M Tris-acetate and 50 mM EDTA.

TAE Buffer is the most commonly used buffer for agarose DNA electrophoresis.

Quality Control:

Nuclease-free.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. Przyrodników 3, 80-297 Gdańsk, Poland

NIP 957-07-05-191, KRS 0000202039, www.eurx.com.pl, orders@eurx.com.pl

tel. +48 58 524 06 97, fax +48 58 341 74 23