



## 5 x TBE

*Molecular Biology Grade*

Cat. No.	size
E0230-01	1 l
E0230-02	5 x 1 l

### **Storage Conditions:**

Store at room temperature.

### **Description:**

- 5 x TBE Buffer is a sterile-filtered solution of 445 mM Tris, 445 mM boric acid and 10 mM EDTA.
- TBE Buffer is the most commonly used buffer for DNA and RNA polyacrylamide gel electrophoresis.

### **Quality Control:**

Each lot is assayed for nuclease activity.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

**EURx Ltd. Przyrodników 3, 80-297 Gdańsk, Poland**  
NIP 957-07-05-191, KRS 0000202039, [www.eurx.com.pl](http://www.eurx.com.pl), [orders@eurx.com.pl](mailto:orders@eurx.com.pl)  
tel. +48 58 524 06 97, fax +48 58 341 74 23