



# Agarose tablets

(CAS NO: 9012-36-6)

Cat. No.	Size
E0305-01	50 tablets
E0305-02	200 tablets

**Grade:** Electrophoresis grade

**Batch No.:** 876906

**Report date:** 07.09.2021

Agarose is supplied in convenient tablet form, featuring 10 tablets in one blister pack. Agarose tablets might be used as safe and clean equivalent of agarose powder in all applications.

Each agarose tablet consists 0.5 g agarose, molecular biology grade, ready to be dissolved in the appropriate amount buffer of choice. Using agarose tablets eliminates the weighting step. Process of dissolving, casting and running electrophoresis gels is performed analogically as with agarose powder.

*For the best results: please allow your agarose tablet presoak in buffer for 1-2 minutes before heating.*

## CERTIFICATE OF ANALYSIS

TEST	SPECIFICATION	RESULTS
Appearance:	White tablet -0.5 g (± 0.01 g)	pass
Solubility:	clear colorless solution at 1 g plus 100 ml water	pass
Water:	≤10%	9.0%
Gel strength (1% Gel):	≥1200 g/cm <sup>2</sup>	1250 g/cm <sup>2</sup>
EEO (Electroendosmosis):	≤0.13	0.06
Sulfate:	≤0.15%	0.12%
Gelling point (1.5% Gel):	36°C ± 1.5°C	37°C
Melting point (1.5% Gel):	88°C ± 1.5°C	89°C
DNase:	None	pass
RNase:	None	pass
Protease:	None	pass
Endonuclease/ligase inhibitory factor:	None	pass

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039  
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23