



## 10 mM Ultrapure dNTPs Mix

Cat. No.	size
E0503-01	1 ml
E0503-02	5 x 1 ml

**Storage Conditions:** Store at -20°C.

Ready to use deoxyribonucleotides mixture.

### Description:

- Ready to use in PCR.
- Free of contaminants.
- Reduction of the number of pipetting steps and the risk of reaction set up errors.

### Concentration:

The concentration of each deoxyribonucleotide is 10 mM (10 mM dATP, 10 mM dGTP, 10 mM dCTP, 10 mM dTTP).

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

**EURx Ltd. Przyrodników 3, 80-297 Gdańsk, Poland**

NIP 957-07-05-191, KRS 0000202039, [www.eurx.com.pl](http://www.eurx.com.pl), [orders@eurx.com.pl](mailto:orders@eurx.com.pl)

tel. +48 58 524 06 97, fax +48 58 341 74 23