



5 mM Ultrapure dNTPs Mix

Cat. No.	size
E0505-01	200 µl
E0505-02	1 ml

Storage Conditions: Store at -20°C.

Ready to use deoxyribonucleotides mixture.

Description:

- Ready to use in PCR.
- Free of contaminants.
- Reduction of the number of pipetting steps and the risk of reaction set up errors.

Concentration:

The concentration of each deoxyribonucleotide is 5 mM (5 mM dATP, 5 mM dGTP, 5 mM dCTP, 5 mM dTTP).

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. Przyrodników 3, 80-297 Gdańsk, Poland

NIP 957-07-05-191, KRS 0000202039, www.eurx.com.pl, orders@eurx.com.pl

tel. +48 58 524 06 97, fax +48 58 341 74 23