



Uncut Human Placental DNA

Cat. No.	size
E4111-01	0.5 mg
E4111-02	2 mg

Storage Conditions:

Store at -20°C.

Store the working solution at 4°C.

References:

1. Kang SH, Park CH, Jeung HC, Kim KY, Rha SY, Chung HC. *Int J Mol Med.* 2007 Jun; 19(6):901-8.

Storage Buffer:

10 mM Tris-HCl (pH 7.4), 1 mM EDTA (1 x TE Buffer).

Quality Control:

A₂₆₀/A₂₈₀ Ratio: 1.71

Y-Chromosome: Material manufactured from only male placentas.

All preparations are assayed for contaminating endonuclease and nonspecific RNase and single- and double-stranded DNase activities.